

## Esco Airstream® Class II Biological Safety Cabinets



The Esco Airstream® Class II biological safety cabinet is an effective solution in providing operator, product and environmental protection within laboratories and industrial facilities. With the presence of its DC ECM blower, this Class II biosafety cabinet can provide 70% energy savings compared to an AC motor. Certified to EN 12469, the Esco Airstream® Class II biological safety cabinet also has antimicrobial coating on all its external and internal painted surfaces for improved safety.

- Sentinel™ Gold Microprocessor Controller displays all safety information on one screen
- Airflow sensor alerts the user if airflow is sufficient for greater safety
- RS 232 Serial Interface Port enables sending of operational information to Building Management System (BMS)
- Stable, self compensating airflow
- ISOCIDE™ powder coat inhibits microbial growth on external surface prevents contamination and improves safety
- Side capture zones and negative pressure side walls optimise containment
- Low noise
- Manual sash
- Compact and ergonomically designed; raised arm rest helps prevent grille blocking and create a comfortable working posture
- UV lamps and electrical outlets are available as factory fitted options - please contact Customer Services for a tailored quotation



### E-Series

- Tempered glass side wall (UV absorbing tempered glass, 5mm, colourless and transparent)
- Multi-piece autoclavable work tray
- Single blower fan system
- Single exhaust filter



Code	Alt Ref	Model	Nominal Size, ft	Ext Dims, w x d* x h, mm	Int Dims, w x d x h, mm	Price
<b>CAB7500</b>	2010718	AC2-2E8	2	730 x 753 x 1400	610 x 580 x 660	<b>£6248.00</b>
<b>CAB7502</b>	2010658	AC2-3E8	3	1035 x 753 x 1400	915 x 580 x 660	<b>£7350.00</b>
<b>CAB7504</b>	2010732	AC2-4E8	4	1340 x 753 x 1400	1220 x 580 x 660	<b>£7350.00</b>
<b>CAB7506</b>	2010656	AC2-5E8	5	1645 x 753 x 1400	1525 x 580 x 660	<b>£10266.00</b>
<b>CAB7508</b>	2010657	AC2-6E8	6	1950 x 753 x 1400	1830 x 580 x 660	<b>£10266.00</b>

\* Without armrest



### S-Series

- One piece stainless steel side wall (1.5mm 16 gauge stainless steel, type 304, with 4B finish)
- Single piece stainless steel work tray
- Single blower fan system
- Single exhaust filter



Code	Alt Ref	Model	Nominal Size, ft	Ext Dims, w x d* x h, mm	Int Dims, w x d x h, mm	Price
<b>CAB7510</b>	2010767	AC2-2S8	2	730 x 753 x 1400	610 x 580 x 660	<b>£6811.00</b>
<b>CAB7512</b>	2010721	AC2-3S8	3	1035 x 753 x 1400	915 x 580 x 660	<b>£7833.00</b>
<b>CAB7514</b>	2010711	AC2-4S8	4	1340 x 753 x 1400	1220 x 580 x 660	<b>£7833.00</b>
<b>CAB7516</b>	2010725	AC2-5S8	5	1645 x 753 x 1400	1525 x 580 x 660	<b>£10769.00</b>

\* Without armrest