

Eppendorf MiniSpin® and MiniSpin® Plus

eppendorf

Entry level centrifuges that are powerful and easy to use, but small enough so that each workstation can be equipped with a personal centrifuge. Optimal ventilation flow reduces heating and protects temperature sensitive samples, and automatic lid opening at the end of the run to prevent sample warming and to allow easy access to samples.

- Extremely compact footprint
- Acceleration and deceleration time <13s
- User-friendly digital display for time and speed
- Separate short spin key for fast and comfortable centrifugation
- Soft touch lid closure for ergonomic lid locking
- Large LCD display
- MiniSpin® plus provides automatic RPM/RCF conversion, and extended timer to 99 min or continuous centrifugation
- Supplied with F-45-12-11 rotor which has a capacity of 12 x 1.5/2.0mL tubes



Code	Description	Supplied with	Max Speed, RPM	Max RCF, xg	Timer	Dims, w x d x h, mm	Price
E5452000060	MiniSpin®	F-45-12-11 rotor	800 - 13400	12100	15s - 30 min	225 x 230 x 130	£898.00
E5452000069	MiniSpin®	F-45-12-11 IVD rotor	800 - 13400	12100	15s - 30 min	225 x 230 x 130	£942.00
E5453000060	MiniSpin® plus	F-45-12-11 rotor	800 - 14500	14100	15s - 99 min	226 x 239 x 120	£1143.00
E5453000062	MiniSpin® plus	F-45-12-11 IVD rotor	800 - 14500	14100	15s - 99 min	226 x 239 x 120	£1200.00

Accessories

Code	Description	Capacity	Pack	Price
E5452727007	F-55-16-5-PCR rotor	2 x PCR strips	Each	£246.00
E5452720008	F-45-12-11 rotor	12 x 1.5/2.0mL tubes	Each	£220.00
E5425716001	Adapter for all 1.5/2.0mL rotors	1 x 0.5mL tube or BD Microtainer® (max Ø 8mm)	6	£12.95
E5452702000	Stainless steel lid for F-45-12-11 rotor	—	Each	£75.00

SLS Lab Basics 12 Place Microcentrifuge

LAB
BASICS

Compact and powerful high speed centrifuge, ideal for research applications including sedimentation of cells and viruses and isolation of macromolecules such as DNA, RNA, proteins or lipids. It consists of a large LCD display which provides all visual information and a unique cooling design which minimises rotor heating.

- Whisper quiet operation – less than 56dB
- Automatic door lock
- Automatic internal diagnosis
- Fully autoclavable polypropylene rotor
- Run time: 30 sec – 99 min
- Pulse operation for quick spins
- Sound alert function
- Over temperature detector
- Over speed detector



Code	Speed Range, RPM	Max RCF, xg	Dims, w x d x h, mm	Price
SLS4600	500 - 15000	15100	245 x 255 x 140	£1154.00