



## Eppendorf 5702/5702 R/5702 RH Centrifuges

eppendorf

The 5702 family of low speed, general purpose centrifuges was especially designed to meet the requirement of low to mid throughput research and cell culture labs. They are compact in design so that one can fit comfortably on every lab bench, quiet in operation to benefit your work environment, and provide 6 different rotor options to accommodate any type of tube. Choose the 5702 for standard applications, the refrigerated 5702 R for heat sensitive samples, or the heated and refrigerated 5702 RH for centrifugation under biological conditions to improve cell viability for more accurate results in downstream applications.

- Very compact footprint to fit on every lab bench
- Very quiet operation to benefit your environment
- Low access height for simple loading and unloading of samples
- SOFT brake function for slow acceleration and braking
- "At set RPM" function starts the timer when the selected rotational speed is reached for reproducible centrifugation runs
- Key lock to prevent unintentional adjustments
- Stainless steel rotor chamber is rust free and easy to clean
- Electronic imbalance detection for maximum safety
- Timer: 0 – 99 min, with continuous run function
- FastTemp® function for fast pre-cooling (5702 R/5702 RH)
- Continuous cooling maintains a constant temperature after the run has been completed (5702 R/5702 RH)
- ECO shutoff engages after 8 hours of non-use to reduce energy consumption and extend compressor life (5702 R/5702 RH)
- Active heating guarantees high temperature accuracy during the entire centrifugation cycle (5702 RH)

	5702	5702 R	5702 RH
Speed Range, RPM	100 – 4400		
Max RCF, xg	3000		
Cooling	No	Yes	Yes
Heating	No	No	Yes
Temperature Control Range, °C	—	–9 to 40	–9 to 42
Dims, w x d x h, mm	320 x 400 x 240	380 x 580 x 260	
Height with Open Lid, mm	530	600	
Weight w/o Accessories, kg	20	36	
Max Power Consumption, W	200	380	

Code	Description	Price
E5702000065	5702 centrifuge (without rotor)	£1569.00
E5703000060	5702 R refrigerated centrifuge (without rotor)	£3720.00
E5704000060	5702 RH refrigerated/heated centrifuge (without rotor)	£4204.00

### Rotors and Accessories

Code	Description	Capacity	Max RCF, xg	Max Speed, RPM	Pack	Price
E5702700002	A-8-17 rotor	8 x 15mL tubes	2800	4400	Each	£238.00
E5702701009	Rubber mat inserts for A-8-17 rotor	15mL round bottom tubes	—	—	8	£15.09
E5702702005	Adapters for A-8-17 rotor	15mL conical tubes	—	—	8	£65.00
E5702704008	F-35-30-17 rotor with 30 x steel sleeves for 15mL vessels, 20 x conical tube adapters, 30 x rubber mats for round bottom tubes	20 x 15mL conical tubes/30 x 15mL glass tubes	2750	4400	Each	£583.00
E5702705004	F-35-30-17 rotor with 10 x steel sleeves for 15mL vessels, 10 x conical tube adapters, 10 x rubber mats for round bottom tubes	20 x 15mL conical tubes/30 x 15mL glass tubes	2750	4400	Each	£388.00
E5702740004	A-4-38 swing bucket rotor without buckets	—	2750	4400	Each	£481.00
E5702720003	A-4-38 swing bucket rotor with 4 x 100mL round buckets	4 x buckets for 1.5 - 100mL vessels	3000	4400	Each	£882.00
E5702722006	Round buckets for A-4-38 rotor	100mL	—	—	4	£497.00
E5702709000	Buckets for A-4-38 rotor	90mL	—	—	4	£300.00