

Eppendorf 5910 Ri Refrigerated Centrifuges Continued

eppendorf

Code	Description	Capacity	Max RCF, xg	Max Speed, RPM	Pack	Price
E5910702000	Adapter for S-4x400 rotor	7 x Eppendorf Tubes® 5mL and 15mL conical	—	—	2	£207.00
E5910703007	Adapter for S-4x400 rotor	15 x 4-8mL tubes	—	—	2	£207.00
E5910704003	Adapter for S-4x400 rotor	11 x 7.5-12mL tubes	—	—	2	£207.00
E5910705000	Adapter for S-4x400 rotor	1 x 175-250mL bottle	—	—	2	£207.00
E5910706006	Adapter for S-4x400 rotor	1 x 400mL wide-neck bottle	—	—	2	£132.00
E5910708009	Adapter for S-4x400 rotor	26 x 1.5/2mL tubes	—	—	2	£207.00
E5910709005	Adapter for S-4x400 rotor	8 x 9mL tubes	—	—	2	£207.00
E5910711000	Adapter for S-4x400 rotor	17 x 12x75mm tubes	—	—	2	£207.00
E5910712006	Adapter for S-4x400 rotor	3 x tubes 50mL skirted, Centriprep and Oak Ridge 30mL	—	—	2	£207.00
E5910714009	Adapter for S-4x400 rotor	1 x 175-225mL conical bottle	—	—	2	£207.00
E5910700008	Aerosol-tight caps for S-4x400 rotor	—	—	—	2	£161.00
E5910710003	Spare sealings for aerosol-tight caps for S-4x400 rotor	—	—	—	2	£99.00
E5895171001	S-4x500 rotor without buckets	—	3220	4000	Each	£1060.00
E5895170005	S-4x500 rotor with 4 x 500mL rectangular buckets	—	3220	4000	Each	£2138.00
E5895155006	FA-30x2 rotor with aerosol-tight QuickLock® lid	30 x 1.5/2mL tubes	20984	13700	Each	£841.00
E5895160000	F-48x15 rotor with 48 x steel sleeves and adapters	40 x 15mL conical tubes	5005	5500	Each	£841.00



Eppendorf 5920 R Refrigerated Centrifuges

eppendorf

Centrifuge 5920 R delivers extraordinary high capacity in a very compact and ergonomic product design. Featuring a superior capacity of up to 4 x 1000mL, a powerful state of the art refrigeration system and advanced temperature management, the 5920 R is the ideal instrument for high throughput applications. The high capacity buckets of the S-4x1000 rotor provide an exceptional high capacity they also allow both tubes and plates to be used in the same rotor bucket.

- Aerosol tight Eppendorf QuickLock® caps and lids allow fast and ergonomic lid locking
- Centrifuge lid with soft touch lid closure
- Quiet operation to improve your work environment
- Compact footprint saves valuable bench space
- Automatic rotor recognition and imbalance detection for maximum operational safety
- FastTemp® function for quick pre-cooling
- FastTemp pro® function allows for automatic pre-cooling based on pre-programmable time and date
- Continuous cooling maintains a constant temperature after the run has been completed – your samples stay cool
- ECO shut off engages after 8 hours (adjustable) of non use to reduce energy consumption and to extend compressor life
- Dynamic compressor control (DCC) technology for precise cooling performance
- Timer: 10s to 99h 59min, with continuous run function, short spin

Code	Description	Supplied with	Price
E5948000060	5920 R refrigerated centrifuge (without rotor)	—	£9301.00
E5948000260	5920 R refrigerated centrifuge	S-4x1000 rotor, high-capacity buckets, 15mL/50mL conical adapters	£11720.00
E5948000365	5920 R refrigerated centrifuge	S-4x1000 rotor, round buckets, 15mL/50mL conical adapters	£11320.00
E5948000460	5920 R refrigerated centrifuge	S-4x1000 rotor, plate/tube buckets, 15mL/50mL conical adapters	£11320.00
E5948000560	5920 R refrigerated centrifuge	S-4x750 rotor, round buckets, 15mL/50mL conical adapters	£11520.00
E5948000961	5920 R refrigerated centrifuge	S-4xUniversal-Large rotor, universal buckets, 5mL/15mL/25mL/50mL conical adapters	£12293.00