



## Eppendorf Cell Culture and Clinical Centrifuge Packages

eppendorf

The Eppendorf centrifuge packages offer a selection of centrifuges supplied with rotors and adapters to meet specific needs.

**Cell Culture Packages:** optimised for swing bucket rotor applications in standard 15mL and 50mL Falcon® tubes and the new swing bucket rotors increase the capacity of Centrifuges 5804/5804 R to 4 x 250mL and Centrifuges 5810/5810 R to 4 x 750mL. They include:

- Centrifuge 5702, 5804/ 5804 R or 5810/ 5810 R
- Swing bucket rotor including buckets
- 4 x adapters for 15mL Falcon® tubes (2 sets of 2)
- 4 x adapters for 50mL Falcon® tubes (2 sets of 2)

**Clinical Packages:** optimised for Centrifuge 5702 R with adapters for 13mm and 16mm blood collection tubes or urine tubes. They include:

- Centrifuge 5702 R
- Swing bucket rotor A-4-38 including round buckets
- 4 x adapters for 13mm blood collection tubes or urine tubes (2 sets of 2)
- 4 x adapters for 16mm blood collection tubes or urine tubes (2 sets of 2)

Cell Culture Centrifuge Packages	Number of Vessels per Adapter/Rotor 15mL Conical	Number of Vessels per Adapter/Rotor 50mL Conical
5702/5702R	2/8	1/4
5804/5804R + rotor A-4-44	4/16	1/4
5804/5804R + rotor S-4-72	8/32	4/16
5810/5810R + rotor A-4-62	9/36	4/16
5810/5810R + rotor A-4-81	12/48	5/20
5810/5810R + rotor S-4-104	14/56	7/28
5810/5810R + rotor S-4-104	23/92 13mm tubes	20/80 16mm tubes

Code	Description	Supplied with	Price
E5702000365	5702 centrifuge	A-4-38 swing-bucket rotor, round buckets, 15mL/50mL conical adapters	£2579.00
E5703000360	5702 R refrigerated centrifuge	A-4-38 swing-bucket rotor, round buckets, adapters for 13mm and 16mm blood collection or urine tubes	£4629.00
E5804000360	5804 centrifuge	A-4-44 swing-bucket rotor, rectangular buckets, 15mL/50mL conical adapters	£4420.00
E5804000565	5804 centrifuge	S-4-72 rotor with adapters for 15mL/50mL conical tubes	£5312.00
E5805000365	5804 R refrigerated centrifuge	A-4-44 swing-bucket rotor, rectangular buckets, 15mL/50mL conical adapters	£6734.00
E5805000665	5804 R refrigerated centrifuge	S-4-72 rotor, adapters for 15mL/50mL conical tubes	£7400.00
E5810000365	5810 centrifuge	A-4-62 swing-bucket rotor, rectangular buckets, 15mL/50mL conical adapters	£5790.00
E5810000460	5810 centrifuge	A-4-81 swing-bucket rotor, rectangular buckets, 15mL/50mL conical	£6510.00
E5810000560	5810 centrifuge	S-4-104 rotor, adapters for 15mL/50mL conical tubes	£6972.00
E5811000360	5810 R refrigerated centrifuge	A-4-62 swing-bucket rotor, rectangular buckets, 15mL/50mL conical adapters	£8433.00
E5811000465	5810 R refrigerated centrifuge	A-4-81 swing-bucket rotor, rectangular buckets, 15mL/50mL conical adapters	£9152.00
E5811000660	5810 R refrigerated centrifuge	S-4-04 rotor, adapters for 15mL/50mL conical tubes	£9614.00
E5810000665	5810 centrifuge	S-4-104 swing-bucket rotor, adapters for 13-16mm round-bottom tubes	£6971.00
E5811000865	5810 R refrigerated centrifuge	S-4-104 swing-bucket rotor, adapters for 13/16mm round bottom tubes	£9614.00