

Thermo Scientific™ SL 8 Series Centrifuges



Perform high capacity work in a compact design. Available in both ventilated and refrigerated models, featuring innovative technologies with a smart, simple interface and the flexibility to meet your clinical and research processes and protocols.

- Compact footprint
- Secure, push button Auto-Lock™ rotor exchange in as little as 3 seconds delivers trouble free rotor installation and removal and clear chamber access for cleaning convenience
- Biocontainment sealing options, including ClickSeal™ lids for glove friendly, one handed operation
- One-touch operation with pre-saved protocols and password protection
- Optional indicators at end of run, including automatic lid opening, flashing display or audible signal
- Ensure a peaceful laboratory with low noise levels
- One-touch operation with pre-saved protocols, glove and detergent friendly interface
- Highly visible backlit display for easy reading of parameters across the lab
- Multilingual instructions
- 4 protocols via direct access key, with naming option and password protection available



	SL 8	SL 8 R
Cooling	No	Yes
Temperature Control Range, °C	—	-10 to 40
Max Speed, RPM	16000	17850
Max RCF, g	24328	30279
Dims, w x d x h, mm	370 x 480 x 310	460 x 670 x 320
Max Power Consumption, W	310	750

Code	Description	Price
75007221	SL 8 centrifuge (without rotor)	£2481.00
75007224	SL 8 R refrigerated centrifuge (without rotor)	£5842.00

Rotors

Code	Description	Capacity	Max RCF, xg	Max Speed, RPM	Price
75005600	MT-12 swinging bucket rotor, 90°, R _{max} 87mm	12 x 2mL tubes	6438	13000	£2335.00
75005701	TX-150 swinging bucket rotor without buckets, 90°, R _{max} 144mm	—	3260	4500	£1021.00
75005706	M-10 microplate rotor without buckets, 90°, R _{max} 119mm	—	2576	4400	£1333.00
75005709	HIGHConic III rotor, 45°, R _{max} 120mm	6 x 50mL tubes	10155 (ventilated)/ 12108 (refrigerated)	8700 (ventilated)/ 9500 (refrigerated)	£2001.00
75005715	MicroClick 24 x 2 rotor, 45°, R _{max} 85mm	24 x 1.5/2mL tubes	24328 (ventilated)/ 30279 (refrigerated)	16000 (ventilated)/ 17850 (refrigerated)	£1252.00
75005719	MicroClick 30 x 2 rotor, 45°, R _{max} 99mm	30 x 1.5/2mL tubes	21694	14000	£1201.00
75005720	8 x 8 PCR strip rotor, 45°, R _{max} 71mm	8 x 8 PCR strips	17860	15000	£1163.00
75005733	Hematocrit rotor with linear scale, 90°, R _{max} 85mm	24 capillaries	16810	13300	£950.00

Rotor Accessories

Code	Description	For use with	Pack	Price
75005731	Replacement lid	HIGHConic III rotor	Each	£212.00
75003058	Replacement sealing O-rings	HIGHConic III rotor	Each	£60.00
75005802	Adapter for 38mL round bottom tubes	HIGHConic III rotor	2	£124.00
75005803	Adapter for 16mL round bottom tubes	HIGHConic III rotor	2	£124.00
75005804	Adapter for 12mL round bottom tubes	HIGHConic III rotor	2	£124.00
75005805	Adapter for 6.5mL round bottom tubes	HIGHConic III rotor	2	£124.00
75005770	Adapter for 5mL round bottom tubes	HIGHConic III rotor	2	£99.00
75005806	Adapter for 3.5mL round bottom tubes	HIGHConic III rotor	2	£124.00
75005807	Adapter for 1.5/2mL conical microtubes	HIGHConic III rotor	2	£124.00
75005702	Round buckets for 4 x 145mL bottles	TX-150 rotor	4	£636.00
75005703	Conical buckets for 8 x 50mL tubes	TX-150 rotor	4	£900.00