

Thermo Scientific™ SL 8 Series Centrifuges Continued

Thermo
SCIENTIFIC

Code	Description	For use with	Pack	Price
75005707	ClickSeal biocontainment lids	Round buckets in TX-150 rotor	4	£214.00
75005724	Replacement O-rings for lids	ClickSeal biocontainment lids for round buckets in TX-150 rotor	2	£96.00
75005735	Adapter for 100mL round bottom open top tube	Round buckets in TX-150 rotor	4	£312.00
75005736	Adapter for 50mL conical skirted or unskirted tube	Round buckets in TX-150 rotor	4	£312.00
75005737	Adapter for 15mL conical tube	Round buckets in TX-150 rotor	4	£312.00
75005739	Adapter for 5/7mL vac blood collection tube	Round buckets in TX-150 rotor	4	£312.00
75005740	Adapter for 3/5mL tubes	Round buckets in TX-150 rotor	4	£312.00
75005743	Adapter for 1.5/2mL tubes	Round buckets in TX-150 rotor	4	£312.00
75005744	Adapter for 25mL universal/flat tubes	Round buckets in TX-150 rotor	4	£312.00
75005808	Adapter for 15mL conical tubes	Conical buckets in TX-150 rotor	2	£230.00
75005723	Unsealed buckets for 4 x standard or 2 x deepwell plates	M10 microplate rotor	2	£747.00
75005725	ClickSeal biocontainment lids	MicroClick 24 rotor	Each	£227.00
75005730	ClickSeal biocontainment lids	8 x 8 PCR strip rotor	Each	£212.00
75005752	Adapter for 0.2mL PCR tubes	Microliter 48 x 2, MicroClick 30 x 2, MicroClick 24 x 2 rotors	30	£225.00
75005753	Adapter for 0.5mL microtubes	Microliter 48 x 2, MicroClick 30 x 2, MicroClick 24 x 2 rotors	30	£239.00
75005754	Adapter for 0.25mL microtubes	Microliter 48 x 2, MicroClick 30 x 2, MicroClick 24 x 2 rotors	30	£248.00

Thermo Scientific™ SL 1 Plus Series Centrifuges

Thermo
SCIENTIFIC

Ergonomic design, standard high contrast user interface for up to 6 saved programs with simple push button operation for routine applications. With capacity up to 1.6L, including 76 x 5/7mL blood tubes and 16 x 50mL conical tubes and Auto-Lock for fast rotor exchange.

- Intuitive controls simplify operation and provide detailed information on a wide range of operating parameters and processing status for optimal performance
- Proprietary motorised lid latch enables one finger downward motion for easiest closing and locking of centrifuge lid, even on tall benches
- Auto-Lock™ secure locking system allows easy push-button installation and exchange of rotors, saves time switching between applications, and allows easy access to rotor chamber for quick cleaning, a healthier working environment and longer unit life
- ClickSeal™ bucket sealing system
- SMARTSpin™ technology for advanced rotor management system maximises acceleration, braking and residual load imbalance for each rotor and bucket option
- Compact dimensions optimise work area and ergonomic height and curved design ensure easy sample loading and unloading and easy cleaning
- Up to 40% energy consumption savings on industry standard protocols, such as standard blood separation or conical tubes processing



	SL 1 Plus	SL 1R Plus
Cooling	No	Yes
Temperature Control Range, °C	—	–10 to 40
Max Speed, RPM		15200
Max Speed, xg		25830
Dims, w x d x h, mm	566 x 660 x 362	625 x 660 x 362
Weight w/o Rotor, kg	57.5	91.5

Code	Description	Price
75009600	SL 1 Plus centrifuge (without rotor)	£4187.00
75009630	SL 1R Plus refrigerated centrifuge (without rotor)	£7155.00