



A

## SLS Select Indicating Tapes



Autoclavable tape indicates when steam sterilisation has occurred by means of a visible black strip which shows after autoclaving for 3 minutes at 134°C, or 15 minutes at 121°C.

Code	Dims, mm x m	Price
<b>SLS1610</b>	19 x 50	£9.69
<b>SLS1612</b>	25 x 50	£12.55



A

## 3M™ Comply™ Autoclave Indicator Tapes

**3M** Science.  
Applied to Life.™

An exposure monitor that securely seal packs and enables personnel to know at a glance that the packs have been exposed to the sterilisation process. The chemical indicator ink, printed in a diagonal striped pattern, changes colour or turns darker when exposed to the sterilisation process.

- Trusted, high quality adhesive
- Lead free ink means no hazardous waste issues
- Stretchable backing that minimises tape "pop-off" during sterilisation
- Can be written on or labelled with pre-printed labels

Code	Alt Ref	Dims, mm x m	Price
<b>AUT1610</b>	1322-18	18 x 55	£10.17
<b>AUT1612</b>	1322-24	24 x 55	£16.28



A

## 3M™ Thermalog™ Integrating Indicator Strips

**3M** Science.  
Applied to Life.™

Class 5 steam chemical integrating indicators offering immediate, accurate and easy to read method of monitoring the sterilisation process conditions inside each pack. The easy-to-read moving front style provides "Safe" or "Unsafe" result, that has been shown to have a higher reading accuracy than colour interpretation style indicators.

- Can be used in all 121 - 135°C steam sterilisation cycles
- Correlates with a Biological Indicator at three time/temperature relationships
- Meet the performance requirements of ANSI/AAMI/ISO 11140-1:2005 Class 5
- Identifies non-sterile instruments before they enter the sterile field
- Helps identify sterilisation process failures by isolating individual pack problems

Code	Alt Ref	Dims, w x l, mm	Pack	Price
<b>AUT1628</b>	2134MM	101 x 19	250 strips	£150.00

## Just a click away...

Our website aims to improve our customer's online experience. With over 500,000 products available and easy navigation, visit today for more information.

[scientificlabs.co.uk/safety](http://scientificlabs.co.uk/safety) (UK & Northern Ireland)

[scientificlabs.ie/safety](http://scientificlabs.ie/safety) (Republic of Ireland)

[scientificlabs.com/safety](http://scientificlabs.com/safety) (International)

