



## Hettich ROTINA 380/380 R Centrifuges



High performance unit supporting a wide array of accessories for multiple applications in just one rotor. Its size and performance makes it the ideal bench top centrifuge for most clinical, research and industrial laboratories.

- Ideal for medium to HTP labs processing large quantities of blood collection tubes, in swing- out rotors and angle rotors
- Metal housing and lid
- Lid with viewing port, powered one-hand lock, dropping protection and emergency lid lock release
- Automatic rotor recognition
- Imbalance switch-off
- IvD-conform according to directive 98/79/EC
- Max. noise level of 54 dB(A) with rotor 1760
- Easy operation with keypad and control knob
- 98 program memories for more individuality
- 9 individual acceleration and 10 deceleration stages

|                               | ROTINA 380      | ROTINA 380 R    |
|-------------------------------|-----------------|-----------------|
| Cooling                       | No              | Yes             |
| Temperature Control Range, °C | —               | -20 to +40      |
| Max Speed, RPM                | 15000           | 15000           |
| Max RCF, g                    | 24400           | 24400           |
| Dims, w x d x h, mm           | 457 x 600 x 418 | 457 x 750 x 418 |
| Weight, kg                    | 51              | 81              |

| Code           | Alt Ref  | Description                                        | Price           |
|----------------|----------|----------------------------------------------------|-----------------|
| <b>CEN5010</b> | GB-C1701 | ROTINA 380 centrifuge (without rotor)              | <b>£5285.00</b> |
| <b>CEN5012</b> | GB-C1706 | ROTINA 380 refrigerated centrifuge (without rotor) | <b>£8591.00</b> |

## Rotors and Accessories

| Code           | Alt Ref    | Description                                                    | Capacity                                                                                | Max Speed, RPM | Max RCF, xg | Price           |
|----------------|------------|----------------------------------------------------------------|-----------------------------------------------------------------------------------------|----------------|-------------|-----------------|
| <b>CEN5014</b> | GB-A1754   | Swing out 4 place rotor, without suspensions, up top angle 90° | —                                                                                       | 5000           | 4863        | <b>£1333.00</b> |
| <b>CEN5016</b> | GB-A1752   | Bucket for GB-A1754 rotor                                      | 1 x 290mL                                                                               | —              | —           | <b>£205.00</b>  |
| <b>CEN5018</b> | GB-A1751   | Bioseal lid for GB-A1754 rotor bucket                          | —                                                                                       | —              | —           | <b>£62.00</b>   |
| <b>A1783</b>   | GB-A1783-A | Adapter for GB-A1754 rotor                                     | 16 x 1.6 - 7mL blood tubes                                                              | —              | —           | <b>£116.00</b>  |
| <b>GBA1763</b> | GB-A1763-A | Adapter for GB-A1754 rotor                                     | 13 x 4 - 15mL blood tubes                                                               | —              | —           | <b>£116.00</b>  |
| <b>A1771</b>   | GB-A1771-A | Adapter for GB-A1754 rotor                                     | 9 x 15mL conical (Falcon type) tubes                                                    | —              | —           | <b>£116.00</b>  |
| <b>A1772</b>   | GB-A1772-A | Adapter for GB-A1754 rotor                                     | 4 x 50mL conical (Falcon type) tubes                                                    | —              | —           | <b>£116.00</b>  |
| <b>A1779</b>   | GB-A1779   | Adapter for GB-A1754 rotor                                     | 5 x 25/30mL conical (Falcon type) tubes                                                 | —              | —           | <b>£148.00</b>  |
| <b>CEN5020</b> | GB-A1760   | Swing out 2 place rotor, without suspensions, up top angle 90° | —                                                                                       | 4000           | 2397        | <b>£1003.00</b> |
| <b>CEN5022</b> | GB-A1753-A | Bucket for GB-A1760 rotor                                      | 10 x microtitre, 6 x culture, 2 x deep well, and 2 x filter plates, 2 x micronic system | —              | —           | <b>£371.00</b>  |