


Kinesis Vials and Caps Continued



Accessories



Code	Alt Ref	Description	Colour	Material	Price
CHR2186	VR12-50PP	Vial rack 12mm x 50 holes	Blue	Polypropylene	£22.00

Thermo Scientific™ Chromacol Vials and Caps



Snap Cap Vials, 11mm

- High quality 11mm crimp top vials made of superior quality 33 expansion borosilicate clear (Type 1, Class A)
- Dims, h x Ø: 32 x 12mm




Code	Alt Ref	Capacity, mL	Residual Volume, µL	Description	Colour	Writing Area	Pack	Price
VIA3084	03-FIRV	0.3	<3	Fused insert	Clear	Yes	500	£409.00
VIA3086	03-FIRV(A)	0.3	<3	Fused insert	Amber	Yes	500	£452.00
VIA3080	1.5-HRRV	1.5	<4	High recovery	Clear	No	100	£255.00
VIA3076	2-RV	2	<170	—	Clear	Yes	500	£108.00
VIA3078	2-RV(A)	2	<170	—	Amber	Yes	500	£143.00

Snap Cap Closures, 11mm

- Eliminate the need for crimping or decapping tools with 11mm Snap-It caps
- Enlarged open area allows for needle penetration across the entire vial opening
- Easy to apply and easy to remove

Supplied in packs of 500



Code	Alt Ref	Description	Price
VIA3088	11-PSN(B)	Integral PE seal and thinned penetration area	£15.00
VIA3090	11-PSN(B)-T02	Pre-fitted with a PTFE seal	£80.00
VIA3092	11-PSN(B)-TST1	Pre-fitted with a red PTFE/white silicone/red PTFE	£264.00
VIA3094	11-PSN(B)-8RT1	Pre-fitted with a red rubber/clear PTFE seal	£114.00
VIA3096	11-PSN(B)-ST1	Pre-fitted with a white silicone/red PTFE seal	£156.00
VIA3098	11-PSN(B)-ST1X	Pre-slit white silicone/blue PTFE seal	£360.00