



neodisher® N



- For neutralising carried over alkali residues and removing lime residues or other acid soluble residues during automated reprocessing
- Based on phosphoric acid
- Free of surfactants
- Supplied as liquid concentrate

Code	Pack	Price
420133	5L	£68.00



neodisher® Z



- Neutralising agent
- For neutralising carried over alkali residues during automated reprocessing
- Based on organic acids
- Free of surfactants and phosphates
- Supplied as liquid concentrate

Code	Pack	Price
420233	5L	£58.00



National Diagnostics OptiClear



- Non-toxic, non-flammable, biodegradable solvent
- Intended for the removal of wax, pitch, flux, grease, resin and solvent soluble resists from optical and electronic components of glass, ceramic, or metal
- May be substituted directly for toluene, xylene or chlorinated hydrocarbons, such as trichloroethylene (TCE), in batch type or ultrasonic cleaners with no alteration in protocol
- High solvency power, yet it will not affect metal or metal alloy, and a neutral pH, so it will not stain glass

Code	Pack	Price
NAT1458	455mL	£58.00
NAT1460	3.8L	£146.00
NAT1462	18.9L	£747.00



Dishwasher Salt

Code	Pack	Price
SAL1010	25kg	£33.70