

Kimtech™ G3 Sterile Gloves

KIMTECH™

- Sterile (Sterility Assurance Level 10⁻⁶)
- Powder free
- Textured fingertips
- Beaded cuff
- Hand specific
- Comfortable fit suitable for double donning
- Cleanroom ready packaging
- Suitable for double donning
- Length: 305mm
- Category III
- AQL 1.5
- EN ISO 374-1:2016 Type C K (chemical risks)
- EN 374-4:2014 (resistance to degradation by chemicals)
- EN ISO 374-5:2016 VIRUS (protection from viruses, bacteria and fungi)

Supplied in packs of 200 pairs

Nitrile

- Skin friendly, white nitrile
- No natural rubber latex, silicone or powder, reducing the risks of skin irritation for the wearer

Code	Alt Ref	Size	Price
SAF4800	HC61160	6	£576.00
SAF4802	HC61165	6.5	£576.00
SAF4804	HC61170	7	£576.00
SAF4806	HC61175	7.5	£576.00
SAF4808	HC61180	8	£576.00
SAF4810	HC61185	8.5	£576.00
SAF4812	HC61190	9	£576.00
SAF4814	HC61110	10	£576.00



Latex

- Natural rubber latex for tactile sensitivity and comfort

Code	Alt Ref	Size	Price
SAF4850	HC1360S	6	£547.00
SAF4852	HC1365S	6.5	£547.00
SAF4854	HC1370S	7	£547.00
SAF4856	HC1375S	7.5	£547.00
SAF4858	HC1380S	8	£547.00
SAF4860	HC1385S	8.5	£547.00
SAF4862	HC1390S	9	£547.00
SAF4864	HC1310S	10	£547.00



SkinGuard 12 Gloves

- For usage up to ISO 14644-1 Class 5 cleanrooms
- Textured surface for increased grip
- 12" cuff length with beading to reduce ripping
- AQL 1.5
- Easy tear cleanroom packaging
- Non sterile
- Ambidextrous

Code	Size	Price
SAF9396	Small	£35.30
SAF9398	Medium	£35.30
SAF9400	Large	£35.30

