

Thermo Scientific™ Sterilin™ 7mL Bijou Containers

Thermo
SCIENTIFIC

- Ideal for small volume samples, such as paediatric samples and other small-volume clinical applications
- Aseptically manufactured
- Suitable for centrifugation at 7200xg
- Leak tested in accordance with BS EN 14254 Annexe D
- Available pre-filled with boric acid (0.09g) for the preservation of urine samples
- Dims, Ø x h: 22.5 x 50.2mm

Supplied in packs of 700

Code	Alt Ref	Description	Price
CON7520	129BBAC	Boric acid, printed label	£426.00
CON7522	129A	No label	£187.00
CON7524	129B	Plain label	£210.00



Thermo Scientific™ Sterilin™ 25mL Certified Containers

Thermo
SCIENTIFIC

- Clear robust polypropylene (PP)
- Convenient, self standing, conical base
- Free from RNase, DNase, human DNA and pyrogens
- Half turn Quickstart cap for ease of handling
- Available plain or with clear blue graduations and large white marking spot
- Maximum RCF: 9500xg
- Lot number on each container for traceability
- Available aseptically manufactured using controlled, automated cleanroom processes or sterile by irradiation
- Excellent leak free performance - tested in accordance with BS EN 14254 Annexe D
- Dims, Ø x h: 31 x 94mm

Supplied in packs of 400

Code	Alt Ref	Description	Sterility	Price
CON7496	30APPRN	Plain	Non sterile	£203.00
CON7498	30BPPRN	Graduated	Non sterile	£210.00
CON7492	30APPRN/IRR	Plain	Sterile	£246.00
CON7494	30BPPRN/IRR	Graduated	Sterile	£253.00



Over 500,000 products available in a click!

- Fully secure e-commerce capability
- Detailed product information
- Leading brands
- Fast and easy ordering
- Visit the [SLS website to find out more](#)
- Promotions and offers

