



Thermo Scientific™ Sterilin™ 30mL Quickstart Containers **Thermo** SCIENTIFIC

- 'Quickstart' cap offers unrivalled, leak free performance
- Excellent temperature and chemical resistance
- 95kPa compliant and CE marked
- Half turn cap ensures improved end user handling
- Available pre-filled with boric acid (0.4g) for the preservation of urine samples
- Supplied in convenient shelf packs of 50
- Suitable for centrifugation at 3800xg
- Lot number on each container aids full traceability
- Aseptically manufactured
- *In vitro* use only
- Containers with spoon have flow-seal style cap
- Dims, Ø x h: 31 x 94mm

Supplied in packs of 400

Code	Alt Ref	Label	Spoon	Material	Price
CON7500	128A	No label	No	Polystyrene	£128.00
CON7502	128B	Printed label	No	Polystyrene	£128.00
CON7504	128C	Plain label	No	Polystyrene	£128.00
CON7506	128BBAC	Boric acid, printed label	No	Polystyrene	£140.00
CON7510	128SB	Printed label	Yes	Polystyrene	£202.00
CON7512	128SC	Plain label	Yes	Polystyrene	£202.00
CON8000	128A/P	Universal, no label	No	Polypropylene	£125.00
CON8002	128B/P	Universal, printed label	No	Polypropylene	£125.00
CON8004	128C/P	Universal, plain label	No	Polypropylene	£125.00
CON8006	128BBAC/P	Boric acid, printed label	No	Polypropylene	£137.00



Thermo Scientific™ Sterilin™ 60-150mL Containers **Thermo** SCIENTIFIC

- Ideal for sample containment
- Aseptically manufactured
- Leak-tested in accordance with BS EN 14254 Annexe D
- 125AM, 125BM and 125CM product certified 95kPa compliant ensuring suitability for sample transport
- Polystyrene (PS)
- Available with either plastic or metal flow seal caps

Supplied in packs of 300

Code	Alt Ref	Capacity, mL	Dims, Ø x h, mm	Label	Closure Type	Pack	Price
CON7530	125AM	60	44.5 x 61	No label	Metal flow seal cap	300	£143.00
CON7532	125BM	60	44.5 x 61	Printed label	Metal flow seal cap	300	£143.00
CON7534	125CM	60	44.5 x 61	Plain label	Metal flow seal cap	300	£143.00
CON7536	125AP	60	44.5 x 61	No label	Plastic wadded cap	300	£143.00
CON7538	125BP	60	44.5 x 61	Printed label	Plastic wadded cap	300	£143.00
CON7540	125CP	60	44.5 x 61	Plain label	Plastic wadded cap	300	£143.00
CON7550	185AM	100	49.5 x 77	No label	Metal flow seal cap	200	£268.00
CON7552	185BM	100	49.5 x 77	Printed label	Metal flow seal cap	200	£267.00
CON7554	185CM	100	49.5 x 77	Plain label	Metal flow seal cap	200	£268.00
CON7556	185AP	100	49.5 x 77	No label	Plastic wadded cap	200	£268.00
CON7558	185BP	100	49.5 x 77	Printed label	Plastic wadded cap	200	£268.00
CON7560	185CP	100	49.5 x 77	Plain label	Plastic wadded cap	200	£268.00
CON7570	165A	150	55 x 108	No label	Metal flow seal cap	120	£196.00
CON7572	165B	150	55 x 108	Printed label	Metal flow seal cap	120	£196.00
CON7574	165C	150	55 x 108	Plain label	Metal flow seal cap	120	£196.00