

Simport® Scientific SecuriTainer™ I Containers



- Designed for collection, transport and storage of specimens
- Shatter resistant which eliminates issues of leakage and evaporation
- Stackable for easy, safe storage and
- Manufactured from virgin polypropylene.
- Magenta polyethylene lids are ribbed for easy opening
- Almost transparent to allow specimens to be viewed without opening.
- Unique tamper evident screw cap using a tab locking mechanism ensures your peace of mind during transport or storage situations where someone might have manipulated the specimen without your prior knowledge
- 95kpa compliant
- Non sterile
- Functional temperature range: -90 to 100°C



Code	Alt Ref	Capacity, mL	Pack	Price
CON2300	C575-20MA	20	500	£146.00
CON2302	C575-40MA	40	500	£180.00
CON2304	C575-60MA	60	500	£170.00
CON2306	C575-90MA	90	400	£171.00
CON2308	C575-120MA	120	300	£148.00

Simport® Scientific SecuriTainer™ III Containers



- Chemically resistant and shatterproof
- Supplied with a leakproof, polyethylene screw cap, particularly important when transporting hazardous material
- Ideal for transport and storage of urine, sputum and most liquids or particulate samples
- Stackable
- Transparent
- Shatter resistant which eliminates issues of leakage and evaporation
- Manufactured from virgin, high clarity polypropylene
- Functional temperature range: -90 to 100°C
- Non sterile



Code	Alt Ref	Capacity, mL	Pack	Price
CON2310	C577-20W	20	500	£131.00
CON2322	C577-40W	40	500	£168.00
CON2324	C577-60W	60	500	£168.00
CON2326	C577-90W	90	400	£170.00
CON2319	C577-120W	120	300	£184.00

Azlon® Screw Cap Containers



- Wide neck with linerless screw caps
- Polypropylene
- Good chemical resistance
- Ideal for storage of samples, powders, etc
- Autoclavable

Code	Alt Ref	Capacity, mL	Dims, Ø x h, mm	Pack	Price
JAR2004	JWP0030P	30	43 x 43	72	£47.10
JAR2006	JWP0060P	60	53 x 47	48	£40.20
JAR2010	JWP0125P	125	70 x 60	36	£40.30
JAR2016	JWP0250P	250	89 x 66	36	£81.00
JAR2022	JWP0500P	500	89 x 97	24	£58.00
JAR2026	JWP0950P	950	120 x 100	24	£151.00
JAR2028	JWP1200P	1200	120 x 132	24	£140.00

