


Silica Crucibles

- 99.8% SiO₂
- Can be used up to 1050°C continuously, and 1300°C for short periods
- Inert to most substances except alkalis and some metallic oxides
- Glazed and translucent
- Low form
- Supplied without lids
- All dimensions and capacities are nominal

Code	Capacity, mL	Dims, h x Ø, mm	Price
CRU3000	5	19 x 29	£31.00
CRU3004	15	25 x 41	£26.00
CRU3006	25	47 x 28	£34.60
CRU3008	50	37 x 57	£43.00
CRU3010	80	45 x 67	£50.00
CRU3012	150	51 x 82	£56.00



Nickel-Electro Stainless Steel Crucibles



- Stainless steel 304 spec 18/8 non-magnetic, bright finish
- Nominal wall thickness: 0.7mm
- Much harder and last longer than porcelain types
- Easy to clean and shatterproof

Code	Capacity, mL	Ø, mm	£Price
ED123	15	30	£10.00
ED124	25	35	£9.00
ED125	30	40	£12.00



Lids

- Flat with a depression for centering on crucibles
- Stainless steel
- Extension on one side that serves as a handle

Code	For use with	Price
ED163	ED123	£64.00
ED164	ED124	£6.00
ED165	ED125	£7.00



Porcelain Crucibles

Low Form

- Lids and crucibles supplied separately

Supplied in packs of 5

Code	Capacity, mL	Dims, h x Ø, mm	Price
CRU20000	6	19 x 30	£16.20
CRU20020	11	22 x 35	£18.20
CRU20040	15	25 x 40	£19.80
CRU20060	25	28 x 45	£24.10
CRU20080	50	38 x 60	£37.00
CRU20100	85	44 x 70	£48.30

