



Medium Form

- Lids and crucibles supplied separately

Supplied in packs of 5

Code	Capacity, mL	Dims, h x Ø, mm	Price
CRU20200	5	20 x 25	£96.00
CRU20240	12	28 x 35	£17.70
CRU20260	20	32 x 40	£19.20
CRU20280	45	40 x 50	£27.00
CRU20300	80	55 x 63	£36.60



Gooch Form


- Supplied with lid

Code	Capacity, mL	Dims, h x Ø, mm	Price
CRU2042	25	40 x 35	£10.00



Lids

Code	For use with	Pack	Price
CRU2116	CRU20000, CRU20200	10	£31.90
CRU2118	CRU20020, CRU20240, CRU2042	10	£34.10
CRU20640	CRU20020, CRU20240	5	£18.60
CRU20660	CRU20040, CRU20260	5	£23.10
CRU20680	CRU20060	5	£25.00
CRU2070	CRU20280	Each	£11.00
CRU20720	CRU20080	5	£39.90
CRU20760	CRU20100	5	£45.20



Thermo Scientific™ Nalgene™ Benchtop Dewars

Thermo
SCIENTIFIC

- Store samples in liquid nitrogen, dry ice solvent mixtures or ice water for the short term
- Double walled sides and CFC free foam urethane insulation minimise heat leakage so samples stay cold longer
- Pour spouts and bottom grip area for easy removal of used liquids
- High density polyethylene base and HDPE cover with polyethylene coated handle for strong, shatterproof performance
- Vented, insulated cover allows for vapour escape but maintains sample temperature
- Can be used from -196 to 100°C

Code	Alt Ref	Capacity, L	Int Dims, Ø x d, mm	Price
DEW0002	4150-1000	1	95 x 194	£234.00
DEW0004	4150-2000	2	121 x 225	£249.00
DEW0006	4150-4000	4	146 x 287	£272.00
DEW0008	4150-9000	10	197 x 394	£501.00