

Dilvac Stainless Steel Dewars, Tall Form

- Fitted in highest quality stainless steel containers, and therefore rust free, anti-corrosive, all hygienic and strong
- All stainless steel lid with self contained insulation, vent and handle
- Toggle clamps secure the lid to the body for safe transportation
- Designed for safe use and transportation of liquid gases
- Designed with practicality and safety in mind



Code	Alt Ref	Capacity, L	Ext Dims, Ø x d, mm	Int Dims, Ø x d, mm	Price
DEW1010	SS/111	1	116 x 228	85 x 185	£225.00
DEW1012	SS/222	2	142 x 268	106 x 205	£262.00
DEW1014	SS/333	4.5	182 x 350	150 x 280	£326.00
DEW1016	SS/700	7	245 x 335	200 x 285	£818.00

Dilvac Blue Enamelled Steel Dewars, Tall Form

- Low cost utility Dewar flasks
- Safe for use with all liquid gases and ideal for cold traps
- High quality borosilicate glass inner flasks, ultra-low evacuated in attractive blue stove-enamelled steel outer container
- Designed with practicality and safety in mind
- Include an insulated lid, fold down carrying handle and cushioned base to prevent benchtop damage



Code	Alt Ref	Capacity, L	Ext Dims, Ø x d, mm	Int Dims, Ø x d, mm	Price
DEW1030	MS/111	1	116 x 228	85 x 185	£203.00
DEW1032	MS/222	2	142 x 268	106 x 205	£239.00
DEW1034	MS/333	4.5	182 x 350	150 x 280	£277.00
DEW1036	MS/200	10	245 x 445	200 x 375	£703.00

Glass Refills

Code	Alt Ref	For use with	Price
DEW1002	470/111	DEW1010, DEW1030	£154.00
DEW1112	470/222	DEW1012, DEW1032	£148.00
DEW1114	470/333	DEW1014, DEW1034	£176.00
DEW1116	470/700	DEW1016	£596.00