


Simport® Scientific Cryovials® Continued



Capinserts™

Polypropylene colour coded inserts which fit precisely into the Simport® Cryovial® caps for colour identification. Etching on top surface for writing purpose and drain hole on top to prevent LN₂ build-up. Assorted bags include white, blue, red, green and yellow inserts.

Code	Alt Ref	Colour	Pack	Price
CRY2258	T312-1	White	500	£16.60
CRY2260	T312-2	Blue	500	£16.60
CRY2262	T312-3	Red	500	£16.60
CRY2264	T312-4	Green	500	£16.60
CRY2266	T312-5	Yellow	500	£16.60
CRY2270	T312-8	Tan	500	£16.60
CRY2272	T312-9	Grey	500	£16.60
CRY2274	T312-10	Lilac	500	£16.60
CRY2276	T312-11	Burnt orange	500	£16.60
CRY2278	T312-13	Violet	500	£16.60
CRY2280	T312-14	Pink	500	£16.60
CRY2268	T312-7	Assortment	5 bags of 100	£16.60

CLEARLine® CryoGen® Tubes/1D Linear Barcode



Innovative screw caps eliminate leakage and contamination. Both internal and external screw caps are co-moulded with a thermoplastic elastomer (TPV) layer and are 95kPa certified to provide a 100% leak proof seal. This proprietary moulding process eliminates the risk of contamination associated with brands that use separate o-rings.

- Human DNA, DNase, RNase, pyrogen, ATP, PCR inhibitors free certified
- CryoGen® Tubes ClearLine® certified as sterile were treated with a validated method to provide a SAL of 10⁻⁶
- CE-IVD marked
- Certified IATA PI 650
- Manufactured in Clean Room Class ISO 7 (UNI EN ISO 14644-1) Class 10000 (US FED STD 209E)
- Working temperature: down to -196°C
- Printed graduations
- Wide white writing surface for specimen identification
- Star shaped top cavity for automatic decapping
- Standard code 128 barcode on each vial with human readable characters, unique and consecutive
- Innovative screw caps eliminate leakage and contamination
- Medical grade raw materials
- Autoclavable at 121°C steam for 20 minutes
- Centrifuged at: 14000RPM/19000xg



Internal Thread

Supplied in packs of 500 (inner sealed zip bags of 50 with sterilisation indicator)

Code	Alt Ref	Volume, mL	Bottom	Price
CRY1056	CL1ARBIPSTS	1.2	Self standing	£174.00
CRY1058	CL2ARBIPSTS	2	Self standing	£187.00
CRY1060	CL2ARIPSTS	2	Round bottom	£181.00
CRY1062	CL4ARBIPSTS	4	Self standing	£163.00
CRY1064	CL4ARIPSTS	4	Round bottom	£163.00
CRY1066	CL5ARBIPSTS	5	Self standing	£171.00
CRY1068	CL5ARIPSTS	5	Round bottom	£171.00