



A

Heathrow Scientific Microtube Storage Boxes



- Made from polypropylene
- Durable three point hinges and a clasp that locks for secure lid closure
- Featuring moulded grid lines on the lid
- Imprinted coordinates on the base, and imprinted and raised coordinates on the box
- User friendly and stackable
- Black, light tight version available for light sensitive samples
- Autoclavable
- Suitable for storage down to -80°C

Supplied in packs of 5

Code	Alt Ref	Capacity	Colour	Dims, w x d x h, mm	Price
RAC0470	HS120032	50 x 1.5 - 2.0mL tubes	Natural	141 x 92 x 56	£45.00
RAC0472	HS120033	50 x 1.5 - 2.0mL tubes	Assorted (blue, green, purple, yellow, orange)	141 x 92 x 56	£45.00
RAC0484	HS120043	100 x 1.5 - 2.0mL tubes	Black	141 x 151 x 57	£55.00
RAC0486	HS120201	100 x 1.5 - 2.0mL tubes	Blue	141 x 151 x 57	£55.00
RAC0487	HS120202	100 x 1.5 - 2.0mL tubes	Green	141 x 151 x 57	£55.00
RAC0488	HS120203	100 x 1.5 - 2.0mL tubes	Yellow	141 x 151 x 57	£55.00
RAC0489	HS120204	100 x 1.5 - 2.0mL tubes	Orange	141 x 151 x 57	£55.00
RAC0480	HS120042	100 x 1.5 - 2.0mL tubes	Natural	141 x 151 x 57	£55.00
RAC0482	HS120044	100 x 1.5 - 2.0mL tubes	Assorted (blue, green, purple, yellow, orange)	141 x 151 x 57	£55.00



Bravo Storage/Cryo Boxes

- Water repellent fibreboard boxes in a choice of colours
- Two heights: 75 or 50mm
- Available with plastic coating for additional waterproofing
- Does not include dividers
- Dimensions, w x d: 133 x 133mm



Code	Colour	Height, mm	Description	Pack	Price
FRE2130	White	50	Standard	12	£31.80
FRE2132	White	75	Standard	12	£38.90
FRE2246	Red	50	Standard	Each	£6.90
FRE2247	Red	50	Plastic coated	Each	£6.80
FRE2216	Red	75	Standard	Each	£5.90
FRE2217	Red	75	Plastic coated	Each	£7.22
FRE2240	Yellow	50	Standard	Each	£5.33
FRE2241	Yellow	50	Plastic coated	Each	£7.30
FRE2210	Yellow	75	Standard	Each	£5.90
FRE2211	Yellow	75	Plastic coated	Each	£7.20
FRE2242	Blue	50	Standard	Each	£5.50
FRE2243	Blue	50	Plastic coated	Each	£7.50
FRE2212	Blue	75	Standard	Each	£5.90
FRE2213	Blue	75	Plastic coated	Each	£7.80
FRE2244	Green	50	Standard	Each	£5.60
FRE2245	Green	50	Plastic coated	Each	£6.80
FRE2214	Green	75	Standard	Each	£6.50
FRE2215	Green	75	Plastic coated	Each	£7.80