


## CLEARLine SBS Footprint Racks

*ClearLine*

- Stackable
- Engraved transparent lid
- Alphanumeric bottom
- Scanning windows 2D and 2D datamatrix on the bottom
- Linear barcode on rack side
- Designed to be used in the nitrogen gas phase up to -196°C
- Autoclavable at 121°C for 20 minutes
- Tubes will securely lock in each well and will not turn for easy screwing and unscrewing by one hand
- For multiscanner and automatic capper-decapper



A

Code	Alt Ref	Capacity	Dims, w x d x h, mm	Pack	Price
CRY1190	BSM58048L	48 x 2mL tubes	85 x 130 x 55	20	£447.00

## Simport® Scientific Cryostore™ Storage Boxes

*Simport*  
Scientific

- Durable cryogenic storage boxes made from extra strong polycarbonate
- Designed to be used at temperatures between -196 and 121°C
- Autoclavable at 120°C, 15psi (1 bar) for 20 minutes
- Transparent cover allows the user to see the contents of the box
- Cover keyed to the base to prevent misalignment control
- Boxes are colour coded, either with coloured grids on the 25 and 81 place boxes, or coloured bases with the 100 place boxes
- Supplied with innovative vial picker to facilitate easy removal of vials



A

Code	Alt Ref	Capacity	Array	Colour	Dims, w x d x h, mm	Pack	Price
CRY2290	T314-225B	25 x 1 - 2mL	5 x 5	White with blue grid	76 x 76 x 52	8	£54.00
CRY2292	T314-225G	25 x 1 - 2mL	5 x 5	White with green grid	76 x 76 x 52	8	£54.00
CRY2294	T314-225R	25 x 1 - 2mL	5 x 5	White with red grid	76 x 76 x 52	8	£54.00
CRY2296	T314-225Y	25 x 1 - 2mL	5 x 5	White with yellow grid	76 x 76 x 52	8	£54.00
CRY2298	T314-281B	81 x 1 - 2mL	9 x 9	White with blue grid	133 x 133 x 53	4	£62.70
CRY2300	T314-281G	81 x 1 - 2mL	9 x 9	White with green grid	133 x 133 x 53	4	£62.70
CRY2302	T314-281R	81 x 1 - 2mL	9 x 9	White with red grid	133 x 133 x 53	4	£62.70
CRY2304	T314-281Y	81 x 1 - 2mL	9 x 9	White with yellow grid	133 x 133 x 53	4	£62.70
CRY2306	T314-481B	81 x 3 - 4mL	9 x 9	White with blue grid	133 x 133 x 81	6	£91.00
CRY2308	T314-481G	81 x 3 - 4mL	9 x 9	White with green grid	133 x 133 x 81	6	£91.00
CRY2310	T314-481R	81 x 3 - 4mL	9 x 9	White with red grid	133 x 133 x 81	6	£91.00
CRY2312	T314-481Y	81 x 3 - 4mL	9 x 9	White with yellow grid	133 x 133 x 81	6	£91.00
CRY2314	T314-581B	81 x 3 - 5mL	9 x 9	White with blue grid	133 x 133 x 95	5	£77.00
CRY2316	T314-581G	81 x 3 - 5mL	9 x 9	White with green grid	133 x 133 x 95	5	£77.00
CRY2318	T314-581R	81 x 3 - 5mL	9 x 9	White with red grid	133 x 133 x 95	5	£77.00
CRY2320	T314-581Y	81 x 3 - 5mL	9 x 9	White with yellow grid	133 x 133 x 95	5	£77.00
CRY2322	T314-2100B	100 x 1 - 2mL	10 x 10	Blue	133 x 133 x 52	4	£42.40
CRY2324	T314-2100G	100 x 1 - 2mL	10 x 10	Green	133 x 133 x 52	4	£42.40
CRY2326	T314-2100R	100 x 1 - 2mL	10 x 10	Red	133 x 133 x 52	4	£42.40
CRY2328	T314-2100Y	100 x 1 - 2mL	10 x 10	Yellow	133 x 133 x 52	4	£42.40