

SLS Lab Pro Statebourne Biosystem Automated Automatic LN₂ Refrigerators Continued



Archive Series Continued

Code	CRY3325	CRY3332	CRY3334	CRY3330	CRY3336	CRY3338
LN ₂ Liquid Phase Capacity, L	Vapour Unit	Vapour Unit	Vapour Unit	Vapour Unit	Vapour Unit	Vapour Unit
LN ₂ Vapour Phase Capacity, L	81	193	236	260	281	300
Weight (empty), kg	356	592	761	808	872	909
Weight (full), kg	711	1287	1776	1976	2197	2315
Ext. Dims, w x d x h, mm	897 x 897 x 1433	1227 x 1227 x 1440	1608 x 1608 x 1272	1608 x 1608 x 1387	1608 x 1608 x 1497	1608 x 1608 x 1587
Int. Dims, w x d x h, mm	750 x 750 x 785	1067 x 1067 x 727	1406 x 1406 x 615	1406 x 1406 x 730	1406 x 1406 x 840	1406 x 1406 x 903
Vial Capacity, 2mL	Vapour Unit	Vapour Unit	Vapour Unit	Vapour Unit	Vapour Unit	Vapour Unit
Vapour Phase Vial Capacity, 2mL	21000	44200	67100	79300	91500	97600

Code	Alt Ref	Model	Price
CRY3325	9916092	Biosystem Archive 21	£31011.00
CRY3332	9916082	Biosystem Archive 40	£43164.00
CRY3334	9916083	Biosystem Archive 60	£48881.00
CRY3330	9916049	Biosystem Archive 80	£53294.00
CRY3336	9916084	Biosystem Archive 100	£55734.00
CRY3338	9916088	Biosystem Archive 100+	£58102.00

Accessories

Inventory systems available in a variety of configurations, please enquire.

Code	Alt Ref	Description	For use with	Price
CRY3300	5801259	Liquid phase inventory system	Biosystem Access 12	£4710.00
CRY3302	5801260	Vapour phase inventory system	Biosystem Access 12	£4375.00
CRY3304	5801252	Liquid phase inventory system	Biosystem Access 24	£7422.00
CRY3306	5801253	Vapour phase inventory system	Biosystem Access 24	£6182.00
CRY3309	5801655	Vapour phase inventory system	Biosystem Access 50	£15990.00
CRY3317	5801623	Vapour phase inventory system	Biosystem Archive 21	£5242.00
CRY3316	5801296	Vapour phase inventory system	Biosystem Archive 80	£17583.00

SLS Lab Pro Statebourne LN₂ Storage Bio Series



Refrigerators with internal canisters for the storage of straws or cryovials on canes. The Bio series offers two ranges of high performance cryogenic refrigerators.

High Capacity Refrigerators

High capacity refrigerators allow the storage of large quantities of biological material while still offering excellent holding times.

- Numbered index canister locations
 - MDD 93/42 EC Compliant
 - High thermal efficiency
 - Extended hold times
 - High capacity
 - Heavy duty and lockable



5 year warranty on vacuum