

## Azlon® Class B Cylinders with Moulded Graduations, PMP



- Calibrated to class B tolerances
- Stable, hexagonal base
- Conforms to ISO 6706
- Transparent and robust PMP
- Good chemical resistance
- Large pouring spout for easy dispensing
- Autoclavable to 121°C

Code	Alt Ref	Capacity, mL	Graduations, mL	Dims, h x Ø, mm	Price
CYL2000	CT0010P	10	0.2	145 x 14	£3.70
CYL2002	CT0025P	25	0.5	170 x 20	£4.92
CYL2004	CT0050P	50	1	200 x 26	£5.53
CYL2006	CT0100P	100	1	250 x 31	£7.50
CYL2008	CT0250P	250	2	315 x 42	£14.10
CYL2010	CT0500P	500	5	360 x 55	£22.90
CYL2012	CT1000P	1000	10	440 x 66	£31.90
CYL2014	CT2000P	2000	20	535 x 83	£66.00




## Azlon® Class B Cylinders with Moulded Graduations, PP



- Moulded graduations
- Stable, anti-suction hexagonal base
- Calibrated to class B tolerance
- Conforms to ISO 6706
- Ultra clear polypropylene
- Good chemical resistance
- Large pouring spout for easy dispensing
- Autoclavable to 121°C

Code	Alt Ref	Capacity, mL	Graduations, mL	Dims, h x Ø, mm	Price
CYL2100	CP0010P	10	0.2	145 x 14	£3.00
CYL2102	CP0025P	25	0.5	170 x 20	£3.04
CYL2104	CP0050P	50	1	200 x 26	£3.21
CYL2106	CP0100P	100	1	250 x 31	£5.90
CYL2108	CP0250P	250	2	315 x 42	£8.70
CYL2110	CP0500P	500	5	360 x 54	£11.50
CYL2112	CP1000P	1000	10	440 x 66	£17.00
CYL2114	CP2000P	2000	20	535 x 83	£36.40



## Azlon® Class B Cylinders with Printed Graduations, PP



- Enhanced, pin sharp, blue printed graduations
- Calibrated to ISO 6706
- Calibrated to class B tolerances
- Ultra clear polypropylene
- Good chemical resistance
- Stable hexagonal base with anti-suction design
- Autoclavable to 121°C

Code	Alt Ref	Capacity, mL	Graduations, mL	Price
CYL2120	CPB0010P	10	0.2	£3.20
CYL2122	CPB0025P	25	0.5	£3.30
CYL2124	CPB0050P	50	1	£3.40
CYL2126	CPB0100P	100	1	£5.60
CYL2128	CPB0250P	250	2	£8.80
CYL2130	CPB0500P	500	5	£11.60
CYL2132	CPB1000P	1000	10	£15.80
CYL2134	CPB2000P	2000	20	£37.50

