


Simport® Scientific CytoSep™ Cytology Funnels Continued



CytoSep™ Metal Funnel Clip

- Holds the cytology funnel disposable sample chamber against the microscope slide
- Compatible with M964 single and double funnels
- Stainless steel

Code	Alt Ref	Price
CYT4008	M964B	£56.00



Simport® Scientific CytoSep™ All Plastic Cytology Funnels



- Easy to use and improve turnaround time while giving more diagnostic information
- Totally disposable to eliminate time consuming decontamination and cleaning
- Safe to use and reduce the risk of exposure to pathogenic samples
- Decrease the possibility of cross contamination
- Do not need metal clips
- Designed to meet any processing requirement
- Capable of producing high quality thin layer slide preparations
- For sample volumes up to 0.5mL
- Cell deposition area is 6mm in diameter (28mm²)
- Filter cards are pre-attached
- Excellent for scanty specimens



Supplied in packs of 40

Code	Alt Ref	Description	Price
CYT4020	M965-10FW	Single funnel with white filter and cap	£87.00
CYT4022	M965-20FW	Double funnel with white filter and cap	£116.00

Simport® Scientific CytoSep™ MEGA Cytology Funnels



- Provides an easy, efficient and cost effective method of producing high quality thin layer slide preparations
- Eliminates the need for stainless steel slide clips
- Dedicated for convenient preparation of larger volume samples
- Quick and easy to remove
- Large rectangular cell deposition area of 22 x 14.75mm (325mm²) provides for up to 12 times the sample volume (6mL) of the single CytoSep™ cytology funnel
- Produces more cost effective thin layer preparations when compared to other thin layer methods
- Strategically placed baffles inhibit cell settling which results in uniform cell deposition and excellent quality slide preparations
- Minimise user exposure to pathogens while reducing the risk of specimen cross contamination.
- Can prepare both air dried and fixed preparations



Code	Alt Ref	Pack	Price
CYT4024	M965-40	40	£116.00

Cytology Sampling Brushes

Designed to provide good quality samples of endocervical cellular material to be used in diagnostic cytological testing. These brushes are designed to be disposable and ensure minimal patient discomfort with maximum sample quality. Also suitable for taking samples of mucous membranes and skin ulcers.

Code	Pack	Price
CYT1050	100	£40.00

