



Reductase Tubes

Manufactured out of soda glass. Dimensions, h x Ø: 150 x 16mm.

Supplied in packs of 100

Code	Capacity, mL	Marks	For	Price
DAI3000	10	1	Methylene blue dye reduction tests	£69.00
DAI3002	2.5/10	2	Antibiotic tests	£95.00
DAI3004	5/10	2	Phosphatase and COD tests	£98.00



SLS Select Peroxide Test Strips



Peroxide test strips, with a measuring range of 0 to 100mg/L H_2O_2 . Colour change white to blue, 0 - 1 - 3 - 10 - 30 - 100mg/L. CE-marked according to the directive for medical products 93/42/EWG.

Code	Strip Dims, w x l, mm	Pack	Price
DAI2110	6 x 95	100 strips	£35.70

SLS Select Dairy Reagents



A larger selection of dairy reagents is available on the SLS website.



Milk Testing Reagents

Code	Description	Pack	Price
DAI1130	0.5% phenolphthalein solution	1L	£55.00
DAI1128	0.5% phenolphthalein solution	2.5L	£66.00
DAI1126	1% phenolphthalein solution	2.5L	£180.00



Phosphatase Test

Code	Description	Pack	Price
DAI1140	4-Nitrophenyl disodium phosphate	0.15g	£4.79
DAI1142	Carbonate bicarbonate buffer (Aschaffenburg and Mullen phosphatase test buffer)	2.5g	£71.00