

SLS Select Dairy Reagents Continued



Cryoscope Standards

Code	Description	Pack	Price
DAI1150	Standard 422 (-0.408°C)	250mL	£17.00
DAI1152	Standard 530 (-0.512°C)	250mL	£17.00
DAI1156	Standard 621 (-0.600°C)	250mL	£17.00
DAI1158	Cryoscope bath liquid	500mL	£19.09



Titration Solutions

Code	Description	Pack	Price
DAI1168	Potassium thiocyanate 0.05N	1L	£55.00
DAI1164	Potassium thiocyanate 0.05N	2.5L	£66.00
DAI1166	Potassium thiocyanate 0.1N	1L	£60.00
DAI1162	Potassium thiocyanate 0.1N	2.5L	£105.00
DAI1170	Silver nitrate 0.01709N	2.5L	£270.00
DAI1178	Silver nitrate 0.02N	2.5L	£358.00
DAI1176	Silver nitrate 0.025N	1L	£266.00
DAI1174	Silver nitrate 0.0282N	2.5L	£169.00
DAI1172	Silver nitrate 0.85N	1L	£166.00



Kartell Vacuum Desiccators

Kartell
LABVACUUM

- Implosion proof
- Will take vacuum of -740mm Hg with loss not exceeding 20mm Hg in 24 hours
- Transparent polycarbonate cover seals onto a neoprene 'O' ring set into the flange of the PP base
- Within the desiccator is a removable pan for holding a drying agent
- Vacuum release tap/non return valve assembly is included, specially grooved stopper in this assembly allows air, when required, to be readmitted very slowly to avoid disturbing the contents

Code	Alt Ref	Capacity, L	Ø, mm	Height, Ext x Int, mm	Price
DES2000	550	2.15	150	190 x 135	£82.00
DES2002	230	4.35	200	230 x 175	£140.00
DES2004	554	9.20	250	300 x 225	£166.00



Polypropylene Desiccator Discs

Supplied in packs of 10

Code	Alt Ref	For use with	Price
DES2010	551	DES2000	£26.00
DES2012	231	DES2002	£37.00
DES2014	553	DES2004	£39.00