



Eppendorf Varispenser® 2/2x Bottle Top Dispensers

eppendorf

The ideal choice for dispensing aliquots of liquid from supply and reagent bottles. Advanced PFA sealing lip technology enables dispensing of nearly all kinds of acids, solutions or bases and prevents crystal building from media residue by wiping off the inner cylinder wall during dispensing. The extra air vent cap with Luer-cone for the assembly of a microfilter prevents possible contamination from the supply bottle. Depending on the chosen size, Varispenser 2 and Varispenser 2x come with 3 or 5 additional thread adaptors.

- Standard GL 45 thread on all sizes for most common bottle threads
- Varispenser 2x is additionally equipped with a recirculation valve for priming the dispensing system without loss of reagent
- Safety valve balls in discharge valve prevents leakage when discharge tube is not mounted
- Screwable air vent cap for easy assembly of drying tube or microfilter
- Oval-shaped dispenser housing for easy transportation
- PFA seal of the sliding piston prevents seizing up
- Fixation slide on interior toothed track for rapid and accurate volume adjustment
- High chemical resistance of all dispensing components
- Telescopic aspirating tubes for variable adjustment to bottle height
- Fully autoclavable without disassembly

| Code | Description | Volume Range, mL | Supplied with | Price |
|--------------------|----------------|------------------|---|----------------|
| E4966000010 | Varispenser 2 | 0.2 - 2 | Adapters GL 25, GL 28/S 28, GL 32, GL 38, S 40 and telescopic intake tube (125 - 240mm) | £254.00 |
| E4966000029 | Varispenser 2 | 0.5 - 5 | Adapters GL 25, GL 28/S 28, GL 32, GL 38, S 40 and telescopic intake tube (125 - 240mm) | £254.00 |
| E4966000037 | Varispenser 2 | 1 - 10 | Adapters GL 25, GL 28/S 28, GL 32, GL 38, S 40 and telescopic intake tube (125 - 240mm) | £254.00 |
| E4966000045 | Varispenser 2 | 2.5 - 25 | Adapters GL 32, GL 38, S 40 and telescopic intake tube (170 - 330mm) | £338.00 |
| E4966000053 | Varispenser 2 | 5 - 50 | Adapters GL 32, GL 38, S 40 and telescopic intake tube (170 - 330mm) | £371.00 |
| E4966000061 | Varispenser 2 | 10 - 100 | Adapters GL 32, GL 38, S 40 and telescopic intake tube (170 - 330mm) | £489.00 |
| E4967000014 | Varispenser 2x | 0.2 - 2 | Adapters GL 25, GL 28/S 28, GL 32, GL 38, S 40 and telescopic intake tube (125 - 240mm) | £300.00 |
| E4967000022 | Varispenser 2x | 0.5 - 5 | Adapters GL 25, GL 28/S 28, GL 32, GL 38, S 40 and telescopic intake tube (125 - 240mm) | £300.00 |
| E4967000030 | Varispenser 2x | 1 - 10 | Adapters GL 25, GL 28/S 28, GL 32, GL 38, S 40 and telescopic intake tube (125 - 240mm) | £300.00 |
| E4967000049 | Varispenser 2x | 2.5 - 25 | Adapters GL 32, GL 38, S 40 and telescopic intake tube (170 - 330mm) | £401.00 |
| E4967000057 | Varispenser 2x | 5 - 50 | Adapters GL 32, GL 38, S 40 and telescopic intake tube (170 - 330mm) | £434.00 |
| E4967000065 | Varispenser 2x | 10 - 100 | Adapters GL 32, GL 38, S 40 and telescopic intake tube (170 - 330mm) | £544.00 |

BRAND Dispensette® S Bottle Top Dispensers



Whether you are dispensing solvents, acids, alkalis or saline solutions - the Dispensette® S bottle-top dispenser makes it easy, safe and efficient. With both the Dispensette® S and Dispensette® S Organic, you will always have the right dispenser at hand for a wide range of media. For efficient volume adjustment, you can choose between digital or analogue adjustment, or fixed volume. With numerous accessory options, the Dispensette® S also offers efficient solutions for special applications.

- Volume is adjusted quickly and set reliably
- Easy to calibrate and adjust in order to comply with ISO 9001 and GLP guidelines
- Autoclavable at 121°C
- Easy to dismantle for cleaning
- Replaceable filling and discharge valve with safety ball
- The valve block can be rotated 360° so that the bottle label always faces the user for safety