

FlyStuff Drosophila Tubes



These quality Drosophila tubes are precision moulded from premium grade polystyrene, K-Resin or polypropylene and conveniently offered in either tray, reload or bulk packs.

- **Polystyrene:** glass-like clarity, excellent visibility, safe, cost-saving alternative to glass tubes, not autoclavable
- **K-Resin:** glass-like clarity, rugged durability, highly flexible, virtually unbreakable, scratch-resistant, ideal for shipping, not autoclavable
- **Polypropylene:** virtually indestructible, good clarity, cost savings, fully autoclavable



Code	Alt Ref	Material	Description	Dims, h x Ø, mm	Packaging	Pack	Price
FLY1006	32-109RL	Polystyrene	Narrow	95 x 25	Reload	5 trays of 100 tubes	£88.00
FLY1190	32-109	Polystyrene	Narrow	95 x 25	Tray pack	5 trays of 100 tubes	£89.00
FLY1310	32-116	Polystyrene	Narrow	95 x 25	Bulk	500 tubes	£66.00
FLY1306	32-110	Polystyrene	Wide	95 x 28.5	Tray pack	5 trays of 100 tubes	£115.00
FLY1312	32-117	Polystyrene	Wide	95 x 28.5	Bulk	500 tubes	£99.00
FLY1308	32-111	K-Resin	Narrow	95 x 25	Tray pack	5 trays of 100 tubes	£176.00
FLY1314	32-118	K-Resin	Narrow	95 x 25	Bulk	500 tubes	£146.00
FLY1102	32-113RL	Polypropylene	Narrow	95 x 25	Reload	5 trays of 100 tubes	£91.00
FLY1088	32-113	Polypropylene	Narrow	95 x 25	Tray pack	5 trays of 100 tubes	£155.00
FLY1318	32-120	Polypropylene	Narrow	95 x 25	Bulk	500 tubes	£104.00
FLY1030	32-121	Polypropylene	Wide	95 x 28.5	Bulk	500 tubes	£72.00

FlyStuff Drosophila Closures



These closures for Drosophila vials and bottles allow for breathability for the Drosophila while helping to prevent mite contamination.

Flugs®

Bonded dense weave cellulose acetate plugs which are said to provide a better barrier to mites than cotton and rayon, and especially paper lids and foam plugs. Additionally, Flugs offer minimal shedding compared to cotton and rayon as well as reduced evaporation, which keeps food fresh longer.

Supplied in packs of 1000

Code	Alt Ref	For use with	Approx Ø, mm	Price
FLY1198	49-100	Stock bottles	37	£217.00
FLY1206	49-101	Wide plastic vials	27	£132.00
FLY1192	49-102	Narrow plastic vials	23	£119.00
FLY1008	49-103	Narrow glass vials	22	£119.00

