



FlyStuff Drosophila Anaesthesia Apparatus Continued

Flowbuddy Flow Regulator

Add convenience and precision flow control to any CO₂ Drosophila workstation using the toggle switch to instantly start and stop fly work with ease. The CO₂ flow to the 2 barbed outputs is adjusted with the Flowbuddy's knurled knob in LPM, with visual confirmation from the built-in flowmeter. Withstands input pressure up to 100psi, meaning the tank gauges can safely be left open. Even when CO₂ pressure increases or decreases, the flow regulator will maintain a consistent unvarying flow. The benchtop version features a sturdy base, non-slip feet, gun hook, and has a small footprint, whereas the wall mounted version is engineered for simple, sturdy installation (hardware included) to free up bench space. The barbed input and outputs are nickel-plated solid brass and fit 3mm ID tubing.



Code	Alt Ref	Description	Supplied with	Price
FLY1050	59-122B	Benchtop flow regulator	—	£498.00
FLY1038	59-122BC	Benchtop flow regulator complete system	Blowgun, Flypad standard, 1829 x 38.1mm (l x Ø) tubing	£782.00
FLY1010	59-122BCU	Benchtop flow regulator complete system	Blowgun, Ultimate Flypad, 1829 x 38.1mm (l x Ø) tubing	£989.00
FLY1236	59-122W	Wall mounted flow regulator	—	£547.00
FLY1238	59-122WC	Wall mounted flow regulator complete system	Blowgun, Flypad standard, 1829 x 38.1mm (l x Ø) tubing	£748.00
FLY1240	59-122WCU	Wall mounted flow regulator complete system	Blowgun, Ultimate Flypad, 1829 x 38.1mm (l x Ø) tubing	£961.00

Regulators

Designed for working with your Drosophila Anesthesia Network. A solid 3mm ID hose barb tube connector fits FlyStuff standard size tubing. Single-stage regulator made of chrome plated brass.

Code	Alt Ref	Description	Price
FLY1060	59-142	Single stage CO ₂ regulator	£390.00



Bubblers

Prevent your flies from becoming dehydrated while also helping in the reduction of static electricity on your Flypad. Bubbler attachment heads come in convenient sizes to fit standard Erlenmeyer flasks. 2 x 914mm coil tubes included. Flask not included.

Code	Alt Ref	For use with	Supplied with	Price
FLY1080	59-181	1000mL flask	No. 9 rubber stopper	£121.00
FLY1242	59-180	500mL flask	No. 7 rubber stopper	£110.00



Accessories

Code	Alt Ref	Description	Price
FLY1142	59-123	T fitting for 3mm tubing	£9.33
FLY1150	59-124C	Tubing, clear	£2.40
FLY1180	59-124	Tubing, green	£66.00
FLY1246	59-125	Coiled tubing with tails, green	£43.70
FLY1162	59-140	Hose-to-hose connector	£4.92
FLY1174	59-141	Quick connect hose-to-hose connector	£12.10
FLY1220	59-178	1/4" to 1/8" tube reducer straight-through connector	£9.48
FLY1230	59-179	3/16" to 1/8" tube reducer straight-through connector	£8.23
FLY1252	59-209	Flow control fitting for blowgun	£45.40
FLY1254	59-210	Flow control fitting for flypad	£56.00
FLY1226	46-100	Porous plastic sheet, 30.5 x 30.5cm	£51.00
FLY1244	59-205	Fly Station brush	£112.00