

Flystuff Drosophila Cages



Population Cages

Featuring cylindrical extruded acrylic for a corner-free interior. The durable custom sewn cloth end caps are form-fitted and held tight by adjustable hose clamps. One end cap features extended netting for easy, safe access. To close the fly cage, simply tie the netting in a knot. The cage sits securely on a durable plastic cradle. Includes acrylic tube, fitted cloth end caps (one with netting for entry), hose clamps and tube cradle.

Code	Alt Ref	Description	Dims, h x Ø, mm	Price
FLY1260	59-116	Small population cage	366 x 229	£620.00
FLY1258	59-104	Large population cage	610 x 305	£1352.00



Embryo Collection Cages

Made of durable acrylic, these convenient embryo collection cages feature an extended end cap (flexible LDPE) which accommodates a Petri dish. This unique easy-to-use feature eliminates the need for tape, while maintaining cage integrity. Stainless steel mesh (97µm) is fused to the top of the cage to allow ample breathability and eliminate condensation. Must be cleaned with cold soapy water.

Code	Alt Ref	Description	For use with	Price
FLY1218	59-105	Mini embryo collection cage	35mm Petri dishes	£84.00
FLY1212	59-100	Small embryo collection cage	60mm Petri dishes	£107.00
FLY1214	59-101	Large embryo collection cage	100mm Petri dishes	£154.00
FLY1172	59-106	Mini replacement end caps	FLY1218	£23.30
FLY1154	59-102	Small replacement end caps	FLY1212	£26.70
FLY1224	59-103	Large replacement end caps	FLY1214	£34.20



Mesh Baskets

Made of durable PVC, these threaded baskets are quite handy to make your own mesh sieves.

Code	Alt Ref	Size	Int Ø, mm	Price
FLY1184	46-101	Small	190	£13.49
FLY1182	46-103	Midi	254	£6.95
FLY1176	46-102	Large	125	£17.07