



SLS Lab Pro Insulated Eco Drying Cabinets



SLS Lab Pro eco drying cabinets are designed with energy efficiency in mind, utilising thermal insulation and a user programmable digital control, for drying laboratory glassware, surgical instruments or media.

- Energy saving thermal insulation to 5 sides
- Low power consumption
- Easy accessible control system
- Digital temperature controller
- Digital 7 day 24hr timer
- Economical to run / lower running costs
- Dual circulating fans and easy roll mounting wheels on floor standing units
- Corrosion free stainless steel internal chamber
- Antimicrobial powder coat paint to external surfaces
- Temperature range: 30 to 80°C

Code	Type	Capacity, L	Airflow	Door Style	Dims, w x d x h, mm	Price
SL5022	Bench	113	Static	Sliding	765 x 450 x 732	£1503.00
SL5024	Bench	227	Static	Sliding	1000 x 550 x 807	£1843.00
SL5026	Floor standing	426	Fan assisted	Hinged	635 x 650 x 1870	£3525.00
SL5028	Floor standing	885	Fan assisted	Hinged	1168 x 650 x 1870	£4750.00



Genlab Drying Cabinets



Suitable for glassware and equipment, these ovens are constructed from powder coated sheet steel body, available in either insulated or non-insulated versions and are fitted with fixed shelf runners and removable plated wire grid shelves as standard. Available insulated for improved efficiency plus energy and cost savings.

- Temperature range: 30 to ambient +50°C
- Easy clean powder coated body
- Toughened glass doors
- Fan circulation on 500L and 1000L models
- Insulated version reduces energy requirements by 30% in comparison to the non-insulated versions
- Easy to use thermostatic control
- Independent overheat thermostat
- Lockable castors (IWC425 and 885 models)
- All units have a powder coated mild steel exterior and the insulated cabinets have 304 stainless steel interior
- Fitted with 3 shelves

Code	Alt Ref	Description	Capacity, L	Ext Dims, w x d x h, mm	Int Dims, w x d x h, mm	Shelves/ Positions	Weight, kg	Price
DRY1000	DC125	Non insulated	125	740 x 420 x 600	720 x 400 x 430	3/5	40	£879.00
DRY1002	DC250	Non insulated	250	1000 x 500 x 660	980 x 480 x 520	3/6	55	£1274.00
DRY1006	DC500	Non insulated	500	600 x 650 x 1680	580 x 620 x 1380	3/12	75	£3012.00
DRY1014	DC1000	Non insulated	1000	1180 x 650 x 1680	1160 x 620 x 1380	3/12	130	£3939.00
DRY1020	IWC100	Insulated	100	740 x 420 x 600	670 x 370 x 400	3/5	50	£1201.00
DRY1022	IWC200	Insulated	200	660 x 1000 x 500	930 x 450 x 490	3/6	75	£1707.00
DRY1024	IWC425	Insulated	425	1680 x 600 x 650	530 x 590 x 1350	3/12	140	£4027.00
DRY1026	IWC885	Insulated	885	1680 x 1180 x 650	1120 x 590 x 1350	3/12	210	£5213.00