

## Kartell Drying Rack



The back plate of the 72 place support is made from a single mould of 4mm rigid HIPS. The peg sockets have closed ends eliminating leakage and consequent biohazards. The seamless drip channel collects waste which drains via drain tube provided. 72 easily removable, push in pegs which can accommodate items with neck bore in excess of 15mm are supplied. This rack will not stain or rust. Additional racks can be joined together to make a bank of racks. Complete with mounting kit, draining hose, and fixing template.

Code	Alt Ref	Dims, w x d x h, mm	Price
<b>DRA1020</b>	213	450 x 110 x 630	£161.00




## Azlon® Drying Racks



Wall mounted or free-standing polyamide epoxy coated steel frame with 8 x 110mm and 24 x 80mm pegs for holding beakers, flasks, cylinders etc. Overall dims, w x h: 410 x 300mm

Code	Alt Ref	Description	Price
<b>DRA1000</b>	RWK070	Wall mounted	£72.00
<b>DRA1010</b>	RWK070FS	Free-standing	£105.00



## Astell Nelldorn Drying Rack



Rack for drying flasks, cylinders etc., can be bench or wall mounted. Compact metal frame with base mounted heater and fan which blows hot filtered air through 60 stepped plastic nozzles of varying length. With heat controls and overheat cut-out. For 220/240v, single phase supply, 2.1kW. Variable 120 minute timer. See the SLS website for spare air filters and nozzles.

Code	Alt Ref	Dims, w x d x h, mm	Capacity	Price
<b>DRY5030</b>	DCU012	490 x 190 x 430	60 items	£1679.00



## Just a click away...

Our website aims to improve our customer's online experience. With over 500,000 products available and easy navigation, visit today for more information.

[scientificlabs.co.uk/safety](http://scientificlabs.co.uk/safety) (UK & Northern Ireland)  
[scientificlabs.ie/safety](http://scientificlabs.ie/safety) (Republic of Ireland)  
[scientificlabs.com/safety](http://scientificlabs.com/safety) (International)

