




SLS Select Blender Bags



- Ideal for use in any blender
- Double size (4mm) seal ensures all faults are totally visible, giving the added security of no leakage
- Tested in production at eight times the normal usage
- Gamma irradiated
- Standard bags are made from irradiated virgin food grade polythene
- Eco bags are manufactured from a new polymer which gives the bag the same strength as a standard bag, whilst being over 20% thinner
- Curved bags have a curved lower seal for Circulator® machines
- Strainer bags have a full width filter suspended inside with a mesh size guaranteed to allow the free movement of organisms through it, whilst stopping debris and fat through, making it easier to pipette the sample for testing

Code	Description	For use with	Thickness, µm	Dims, w x l, mm	Pack	Price
SLS3100	Standard	80mL models	65	100 x 150	2000	£79.00
SLS3102	Standard	400mL models	65	180 x 300	500	£72.00
SLS3104	Standard	3500mL models	70	380 x 500	250	£113.00
SLS3103	Curved	400mL circulator® models	65	180 x 300	500	£71.00
SLS3108	Strainer	400mL models	65	180 x 300	200	£112.00
SLS3109	Strainer curved	400mL circulator® models	65	180 x 300	200	£100.00
SLS4102	Eco	400mL models	50	180 x 300	2500	£208.00



Cable Ties

- White PVC

Code	Pack	Price
BA61152	100 ties	£3.12



SLS Lab Pro SM Series Micro Analytical Balances



Developed for high end micro weighing with the automatic elimination of static electricity at a competitive price. A built in fanless ioniser eliminates static and a large weighing chamber with interlockable slide doors allows operation with either hand.

- One touch automatic calibration
- Internal data memory function for up to 200 sets of data
- GLP/GMP/GCP/ISO compliant output
- Built in temperature, humidity, and atmospheric pressure sensors
- RS-232C and Quick USB interfaces as standard
- AD-1688 weighing data logger as standard

Code	Capacity, g	Readability, mg	Pan Dims, mm	Price
SLS5500	22	0.001	Ø 25	£11385.00
SLS5502	5.1/22	0.001/0.01	Ø 25	£10103.00
SLS5504	250	0.01	Ø 90	£5481.00
SLS5506	220	0.1	Ø 90	£3593.00
SLS5508	320	0.1	Ø 90	£3899.00
SLS5510	520	0.1	Ø 90	£4178.00

Accessories

Code	Description	Pack	Price
BAL4651	LAN ethernet interface	Each	£490.00
BAL4656	Tabletop breeze break	Each	£1942.00
BAL4657	Transparent plastic working cover for display and keys	5	£49.00