

Hanna 200 Series Bench Meters



Simple and straightforward bench pH meters designed for every-day use by technical and non-technical staff alike. A two button calibration procedure with pre-programmed buffers, and step by step calibration guide, along with a stability indicator help the user to gain accurate results simply and efficiently.

- Extra wide LCD display
- Calibration expiration reminder helps ensure the meter is in optimum working order
- Sturdy wide-based casing with easy to use buttons for quick efficient testing
- Automatic temperature compensation
- Measurement memory / memory recall function
- Supplied with HI-1131B pH electrode, HI-7662 temperature probe, HI-76404N electrode holder, starter pack of buffer and cleaning solutions, 12 Vdc adapter, and instructions



Specifications	
pH Range	-2.00 to 16.00, ± 0.01
pH Resolution	0.01
Temperature Range, °C	-20.0 to 120.0, ± 0.4 (excluding probe error)
Temperature Resolution, °C	0.1
mV Range (PHM4204 only)	$\pm 399.9 \pm 0.2$; $\pm 2000 \pm 1$
mV Resolution (PHM4204 only)	0.1; 1
Dimensions, mm	235 x 222 x 109
Weight, kg	1.1; kit with holder 2.5

Code	Alt Ref	Description	Price
PHM4202	HI-2210-02	Benchtop pH and temperature meter	£465.00
PHM4204	HI-2211-02	Benchtop pH, temperature and mV meter	£524.00

Accessories

Code	Alt Ref	Description	Price
PHE1150	HI-1131B	Replacement pH electrode	£152.00
PHE1152	FC200B	pH electrode for food use	£200.00
PHM4123	HI-3230B	Platinum ORP probe	£179.00

Sustainable laboratory solutions!

We've planted 2500 trees to offset the carbon footprint of this catalogue by double the amount of CO₂, making it climate positive.

Look out for the following logos on our literature:

