

Sartorius Secura® Micro Balance

SARTORIUS

Secura® gives you the security and peace of mind of knowing that you have done everything right. Besides providing highly accurate weighing results and operating convenience, Secura® also features built-in protection systems for complete reliability and regulatory compliance – the safe and secure way. Real time guidance prompts for leveling, automatic internal adjustment, and 100% traceable, clear documentation with sample and batch identifiers make your lab work more efficient.

- Integrated protection systems for minimising your risk of obtaining incorrect results
- SQmin for automatic detection and active monitoring to ensure you are within the acceptable USP operating range for reliable results
- Fully automatic internal adjustment function, isoCAL, for the highest possible accuracy at all times
- Cal Audit Trail for easy, traceable documentation of all calibration and adjustment data
- Safe and clear documentation by sample and batch identifiers for keeping perfect track of weighing data
- GLP/GMP compliant printouts
- Standard password security features for protection against tampering and unauthorised access
- Integrated PC direct function and plug and play connectivity to a Sartorius printer



Code	Alt Ref	Readability, mg	Capacity, g	Pan Dims, mm	Dims, w x d x h, mm	Price
BAL3872	SECURA26-1S	0.002	21	Ø 50	218 x 376 x 316	£19642.00

Sartorius Quintix® Semi-Micro Balances

SARTORIUS

Quintix® makes laboratory work more enjoyable. The intuitive touch screen in combination with the most precise weighing technology make Quintix® laboratory balances the only choice to achieve the best results. Benefit from improved workflows by easy to use application programs. If you need to generate GLP compliant documents, you will save considerable time and effort with "plug and work" connectivity.

- Intuitive operation of the balance without having to read a manual thanks to the new user interface and colour touch screen
- Integrated application programs and intuitive user guidance for greater efficiency in your daily lab routine
- Highest accuracy and convenience ensured at all times by the isoCAL fully automatic internal adjustment feature
- Easier record keeping using the balances directly integrated PC functions and plug and play connectivity to a Sartorius lab printer
- Pan dims: Ø 80mm



Code	Alt Ref	Capacity, g	Readability, mg	Price
BAL1160	QUINTIX35-1S	30	0.01	£4837.00
BAL1162	QUINTIX65-1S	60	0.01	£5512.00
BAL1164	QUINTIX125D-1S	60/120	0.01/0.1	£5990.00