



Thermo Scientific™ Orion™ ROSS™ pH Electrodes

Thermo
SCIENTIFIC

Choose from a wide variety of Thermo Scientific™ Orion™ ROSS™ glass or epoxy body electrodes, such as: refillable, gel filled, Sure-Flow™, Triode™, semi-micro, micro, flat bulb, spear tipped, and Ultra™ for extended warranty (18 months for gel/24 months for refillable). ROSS™ electrode connectors are also available to fit most titrators.

- Superior measurement stability and fast response
- High accuracy and precision in samples with varying temperatures
- Double junction design means the reference system is isolated from the sample by an intermediate fill solution to make the electrode compatible with samples such as those that contain TRIS, proteins and sulphides
- Reference part of the electrode is maximally protected from the sample and lasts longer than in electrodes with conventional design
- Inert platinum wire does not introduce metal ions into sample and does not deplete over time
- pH range: 0 – 14 (0.01 resolution)

Code	Description	Price
8101BNWP	ROSS® pH half-cell with glass body	£220.00
8102BN	ROSS® combination pH electrode	£277.00
8102BNUWP	ROSS® Ultra® combination pH electrode with extended life	£364.00
8102SC	ROSS® combination pH electrode	£250.00
8103BN	ROSS® combination pH electrode with glass body, semi-micro tip	£344.00
8103BNUWP	ROSS® Ultra® combination pH electrode with glass body, semi-micro tip	£454.00
8103SC	ROSS® combination pH electrode with glass body, semi-micro tip	£507.00
8104BN	ROSS® combination pH electrode with glass body, rugged bulb	£337.00
8104BNUWP	ROSS® Ultra® combination pH electrode with rugged bulb, extended life	£402.00
8104SC	ROSS® combination pH electrode with glass body, rugged bulb	£287.00
8107BNUMD	ROSS® Ultra® pH/ATC Triode with epoxy body, low maintenance gel	£396.00
8107UWMMD	ROSS® Ultra® pH/ATC Triode with epoxy body, low maintenance gel 3M cable	£410.00
8115BNUWP	ROSS® Ultra® combination pH electrode with epoxy body, semi-micro tip	£409.00
8135BN	ROSS® combination pH electrode with epoxy body, flat surface	£347.00
8135BNUWP	ROSS® Ultra® combination pH electrode with epoxy body, flat surface	£422.00
8135SC	ROSS® combination pH electrode with epoxy body, flat surface	£317.00
8156BNUWP	ROSS® Ultra® combination pH electrode with epoxy body	£355.00
8157BNUMD	ROSS® Ultra® pH/ATC Triode with epoxy body	£409.00
8162SC	ROSS® combination pH electrode with 14/15 standard taper	£381.00
8163BNWP	ROSS® combination pH electrode with glass body, spear tip	£445.00
8165BNWP	ROSS® Sure-Flow® combination pH electrode with epoxy body	£383.00
8172BNWP	ROSS® Sure-Flow® combination pH electrode with glass body	£387.00
8302BNUMD	ROSS® Ultra® pH/ATC Triode with glass body	£423.00
8175BNWP	ROSS® Sure-Flow® combination pH electrode with epoxy body, semi-micro tip	£433.00
800500U	ROSS® Sure-Flow® combination pH electrode with epoxy body, semi-micro tip	£240.00
8220BNWP	ROSS® Log R combination pH electrode with glass body, micro	£620.00