

SLS Flowgen multiSUB Midi Horizontal Electrophoresis Systems Continued



Accessories

Code	Description	Price
MS10UV7	UV tray, 100 x 70mm	£66.00
MS10UV10	UV tray, 100 x 100mm	£70.00
MS10UVS	UV gel scoop, 100mm	£45.00
MS10PE	Positive electrode	£142.00
MS10NE	Negative electrode	£142.00
MS10UVDAM	Casting dams	£61.00
MS710FC	multiSUB Mini/Midi flexicaster	£162.00
MS10LG	Adhesive loading guides	£16.00
MS10WP	Viewing platform	£17.00

SLS Flowgen multiSUB Choice Horizontal Electrophoresis Systems



The multiSUB Choice offers a wide degree of versatility. Three tray options are available allowing the choice of one, two or all three gel length options. Maximising comb and tray options allow up to 210 samples to be resolved per gel. The 15cm total run length allows restriction fragment or other close MW sample bands to be easily separated and identified. Speed loading is accomplished using 10, 14, 16, or 18 sample multi channel pipette compatible combs.

- Available with 150 x 70mm, 150 x 100mm, or 150 x 150mm gel trays
- Ideal for restriction fragment analysis, sample prep or checking of high sample numbers
- Run up to 210 samples
- Low buffer volumes
- Multichannel pipette compatible combs for speed loading
- Supplied with 2 x 20 sample combs (1mm), loading guides and dams



Code	Description	Gel Dims, w x l, mm	Max. Capacity	Price
MSCHOICE7	multiSUB Choice 7	150 x 70	70 samples	£385.00
MSCHOICE10	multiSUB Choice 10	150 x 100	140 samples	£470.00
MSCHOICE15	multiSUB Choice 15	150 x 150	210 samples	£384.00
MSCHOICETRIO	multiSUB Choice Trio	150 x 70 & 150 x 100 & 150 x 150	210 samples	£471.00

Combs – 0.75mm Thick

Colour coded: black.

Code	Description	Price
MS151075	Prep 1, marker 1	£44.00
MS152075	Prep 2, marker 2	£44.00
MS154075	Prep 4, marker 2	£44.00
MS1510075	10 sample	£44.00
MS1510MC075	10 sample, multichannel	£44.00
MS1512075	12 sample	£44.00
MS1514MC075	14 sample, multichannel	£44.00
MS1516MC075	16 sample, multichannel	£44.00
MS1518MC075	18 sample, multichannel	£37.80
MS1520075	20 sample	£44.00
MS1530MCSS075	30 sample, multichannel	£34.50
MS1535075	35 sample	£44.00