

## OHAUS Explorer® Semi-Micro Balances



When laboratory work and research requires accuracy to the hundred-thousandth, there is simply no room for error. The Explorer® series of semi-micro balances has been designed with the technology to ensure that your very specific weighing results are accurate.

- Geared for high performance with up to 0.01mg readability and capacities up to 220g
- AutoCal™ Internal Calibration system ensures accuracy of results in applications that require high precision
- Colour touchscreen display that enhances the user experience
- Touchless sensors, automatic draftshield doors (on select models) and ioniser ensure the integrity of results
- SmartText™ 2.0 software, intuitive icon based navigation, 99 item library, 6 user profiles and below minimum sample weight indication – all combine to provide an easy yet efficient weighing experience
- Easy access communication ports including 2 units of USB port and RS232 and an optional fourth port of either RS232 or ethernet, GLP/GMP data output with real time clock, data transfer function
- Pan dims: Ø 80mm



Code	Alt Ref	Model	Capacity, g	Readability, mg	Auto Door	Price
BAL1134	30218998	EX125D	52/120	0.01/0.1	No	£4960.00
BAL1136	30218999	EX125	120	0.01	No	£5255.00
BAL1138	30219000	EX225D	120/220	0.01/0.1	No	£5515.00
BAL1140	30219001	EX225D/AD	120/220	0.01/0.1	Yes	£7325.00
BAL1142	30219002	EX225/AD	220	0.01	Yes	£11710.00

## SLS Lab Pro A-Series Analytical Balances



These top loading analytical balances have a compact footprint. The breeze break has rotary sliding doors that slide behind the balance and can be removed with one touch. These features, combined with the small footprint, make the balances ideal for applications where space is at a premium.

- Compact footprint
- Reverse backlit LCD display for easy to read weighing results
- Statistical calculation function eliminates the need for a stats printer
- Standard RS-232C
- Fast, compact SHS sensor, response time ~2 seconds



Code	Capacity, g	Readability, mg	Calibration	Price
SL55512	102	0.1	External	£908.00
SL55058	102	0.1	Internal	£1103.00
SL55514	152	0.1	External	£1021.00
SL55516	152	0.1	Internal	£1254.00
SL55518	252	0.1	External	£1164.00
SL55060	252	0.1	Internal	£1435.00

### Accessories

Code	Description	Price
BAL4650	Quick USB interface with cable	£98.00
BAL4651	LAN ethernet interface	£490.00
BAL4652	Internal rechargeable NiMH battery pack	£358.00
BAL8020	Transparent plastic working covers for screen and keyboard (5pcs)	£44.00