

2
YEARS

SLS Flowgen powerPRO 3AMP Power Supply



Designed for virtually all high current electrophoresis applications. The PowerPro 3AMP's higher current output capability is perfect for electroblotting units with high-intensity plate electrodes, particularly the omniBLOT maxi, VS20 WAVE and semi dry blotting systems.

- Electro transfers may be performed as timed runs in constant or programmable mode to prevent overheating and buffer depletion
- Shares the built-in protocol and constant parameter features of the other PowerPro models
- Run time extendable to a maximum 999 minutes in constant mode is also useful for overnight transfers undertaken at constant low current in wire electrode systems

Code	Max Voltage, V	Max Current, mA	Max Power, W	Output Pairs	Dims, w x d x h, mm	Price
ELE2604	300	3000	300	5	215 x 335 x 104	£886.00

3
YEARS

Consort Power Supplies

Each model is fully programmable, with up to nine different programs permanently stored at any one time. Each program can have nine parameters, so very complex multistep techniques can be programmed. All units have four parallel, 4mm sockets.

- Method programming
- Clear LCD display
- Automatic crossover
- Data logging facility
- Password protection
- Data acquisition/control software for PC (optional)
- RS232 interface as standard
- Timer 0 - 99.59 hours
- Safety features include ground leaking detection, overload protection, smooth voltage rise, no load detection and isolated communication

Applications	EV2310	EV2650	EV2320	EV3020	EV3610	EV3150	EV3330	EV3620
DNA Sequencing							Y	Y
Flat Bed Isoelectric Focussing						Y	Y	Y
Horizontal Gel	Y	Y	Y	Y	Y	Y	Y	
Long Vertical Gel								Y
Vertical Slab Gel	Y	Y	Y	Y	Y	Y	Y	
Electro-Elution		Y	Y		Y	Y		
Western Blotting				Y				
Semi-Dry Blotting				Y				
Mini Western Blotting	Y							
Mini Semi-Dry Blotting	Y							

Code	Alt Ref	Output voltage, V	Max Current, mA	Max Power, W	Price
ELE9104	EV2230	0 - 1500	0 - 300	150	£1040.00
ELE9106	EV2310	0 - 300	0 - 1000	150	£787.00
ELE9108	EV2320	0 - 3000	0 - 1500	150	£1040.00
ELE9110	EV2650	0 - 600	0 - 500	150	£787.00
ELE9112	EV3020	0 - 300	0 - 2000	300	£1111.00
ELE9114	EV3150	0 - 1200	0 - 500	300	£1111.00
ELE9116	EV3330	0 - 3000	0 - 300	300	£1677.00
ELE9118	EV3610	0 - 600	0 - 1000	300	£1111.00
ELE9120	EV3620	0 - 6000	0 - 150	300	£1677.00