


## Lonza FlashGel™ Electrophoresis Packs

**Lonza**

The FlashGel™ System is the fastest way to separate DNA and the only way to watch DNA migration as it happens. This revolutionary new tool separates DNA in 2 – 7 minutes. Monitor gel runs right at the bench, without UV light. The highly sensitive system allows a 5 – 20x reduction in DNA levels, saving both time and money.

- Get results in 5 minutes or less
- Watch DNA migrate at your bench, in real-time without UV light
- The FlashGel™ system consists of enclosed, disposable, precast agarose gel cassettes and a combination electrophoresis and transilluminator unit
- FlashGel™ cassettes contain precast, prestained agarose gels and buffer – no need for gel preparation, buffer addition or gel staining
- The FlashGel™ dock is an electrophoresis apparatus with a built in transilluminator that provides both separation and detection



Code	Description	Supplied with	Price
<b>LZ57024</b>	RNA starter kit	FlashGel™ cassettes, sample buffer, RNA marker, water (dock sold separately)	£540.00
<b>LZ57026</b>	Starter kit	FlashGel™ dock, pack of 1.2% DNA cassettes, loading dye & marker	£826.00
<b>LZ57062</b>	Device pack	FlashGel™ dock, FlashGel™ camera, FlashGel™ power supply	£1994.00
<b>LZ57065</b>	Camera pack	FlashGel™ dock, FlashGel™ camera	£1590.00
<b>LZ57067</b>	System kit	FlashGel™ dock, FlashGel™ camera, 1.2% DNA cassettes, loading dye & marker	£1937.00
<b>LZ57069</b>	Power supply pack	FlashGel™ dock, FlashGel™ power supply	£1162.00

## Lonza FlashGel™ Cassettes

**Lonza**

Pre-cast gels for use with the Lonza FlashGel™ electrophoresis system.

Supplied in packs of 9

Code	Description	Price
<b>LZ57023</b>	DNA cassettes 1.2%, 12+1 wells	£191.00
<b>LZ57027</b>	RNA cassettes 1.2%, 12+1 wells	£168.00
<b>LZ57028</b>	RNA cassettes 1.2%, 2 x 16+1 wells	£162.00
<b>LZ57029</b>	DNA cassettes 1.2%, 2 x 16+1 wells	£154.00
<b>LZ57031</b>	DNA cassettes 2.2%, 12+1 wells	£158.00
<b>LZ57032</b>	DNA cassettes 2.2%, 2 x 16+1 wells	£158.00
<b>LZ57022</b>	Recovery 2.2%, 8+1 well double tier format	£154.00



## Over 500,000 products available in a click!

- Fully secure e-commerce capability
- Detailed product information
- Leading brands
- Fast and easy ordering
- Visit the [SLS website to find out more](#)
- Promotions and offers

