

## Eppendorf Eporator®

eppendorf

A compact instrument designed for fast and controlled electroporation of bacteria, yeasts and other microorganisms. The instrument features an intuitive operation and user friendly programming of standard methods. Experimental data can easily be exported and documented using its USB port.

- Faster sample handling – simple one button operation
- Maximum safety via integrated cuvette holder
- Compact, space saving design
- Self explanatory operation guided by display indicators
- Fast, freely programmable softkeys for saving and accessing the 2 most frequently used parameter settings
- USB interface for data transfer and GLP supporting documentation



Specifications	
Pulse Voltage, V	200 – 2500
Pulse Form	Exponentially decreasing
Time Constant	Nominal 5ms
Power, W	20
Resistance, $\Omega$	600
Capacitor, $\mu$ F	10

Code	Dims, l x w x h, mm	Weight, kg	Power Supply	Interface	Price
E4309000035	275 x 190 x 125	3.2	100 – 240V, 50 – 60Hz	USB 2.0	£2852.00

## Accessories

Stackable rack with non slip silicone feet and a tilt proof footprint design. Polypropylene for high chemical and physical resistance. Optimised space saving dimensions. Numbered tube positions. Autoclavable.

Code	Description	Pack	Price
E0030119851	30 position cuvette rack	2	£40.50



## SLS Flowgen Electroporation Cuvettes



Designed to maximise molecular electroporation and electrofusion efficiencies for bacteria, yeast, insect, plant and mammalian cells. Universal cuvettes which are compatible with most electroporation systems.

- Consistent uniform pulse generation and improved gene transfer
- Uniform electrode surface and maximises transformation efficiencies
- Cuvette caps engineered to improve aseptic handling techniques and reduce contamination
- Tapered 'V' bottom design reduce sample volumes used
- Available in 1 and 2mm gap sizes with individual colour coded caps
- Sterile and individually wrapped

Supplied in packs of 50

Code	Description	Price
FBR-101	Short, 1mm	£105.00
FBR-102	Short, 2mm	£104.00
FBR-201	Long, 1mm	£107.00
FBR-202	Long, 2mm	£105.00
FBR-204	Long, 4mm	£107.00

