

Buchi R-100 Rotavapors® Continued



Accessories

Code	Alt Ref	Description	Pack	Price
EVA4164	000636	Vacuum seal KD22	Each	£46.00
EVA4170	032339	Vapour duct SJ29/32 with combi-clip (for V and C condensers)	Each	£130.00
EVA4130	048618	Flask holder	Each	£32.80
EVA1095	11059916	Flask holders	5	£116.00



Buchi V-100 Vacuum Pumps



An economical and chemical resistant PTFE diaphragm vacuum pump with a space saving design. It covers the essential needs in vacuum generation. Optimally suited to work with the Rotavapor® R-100.

- Effective vacuum source for limited budgets
- Small footprint saves valuable laboratory space
- Useable for multiple applications
- Lightweight and compact design with integrated carrying handle
- Low maintenance
- Highly chemical resistant thanks to use of quality materials
- Integrated ECO² mode, active after 1hr with suction capacity of 1m³/hr
- All configurations supplied with vacuum hose, power cable, silencer and set hose connection

Specifications	
Suction Capacity, m ³ /h	1.5
Final Vacuum (absolute), mbar	10 (±2)
Max Nominal Speed, RPM	1280
Nominal Speed in ECO ² Mode	70% of nominal speed

Code	Alt Ref	Description	Control	Condensation Assembly	Price
EVA1102	11V100000	V-100	Manual	—	£1572.00
EVA1103	11V100001	V-100	Manual	Secondary condenser with insulation	£2159.00
EVA1104	11V100002	V-100	Manual	Secondary cold trap	£2159.00
EVA1105	11V100010	V-100 with Woulff bottle	Manual	—	£1539.00
EVA1106	11V100011	V-100 with Woulff bottle	Manual	Secondary condenser with insulation	£2327.00
EVA1107	11V100012	V-100 with Woulff bottle	Manual	Secondary cold trap	£2327.00
EVA1098	11V100110	V-100 with Woulff bottle	1-100 interface	—	£2519.00
EVA1099	11V100111	V-100 with Woulff bottle	1-100 interface	Secondary condenser with insulation	£3302.00
EVA1101	11V100112	V-100 with Woulff bottle	1-100 interface	Secondary cold trap	£3302.00

Accessories

Code	Alt Ref	Description	Price
EVA4224	11062291	Silencer for Buchi V-700/V-710 vacuum pumps	£85.00
EVA4223	11062144	Replacement diaphragm for Buchi V-700/V-710 vacuum pumps with fitting tool	£343.00