

Buchi R-300 Rotavapors® Continued

BUCHI

Accessories

Code	Alt Ref	Description	Price
EVA1060	11062186	Vapour duct for SJ 29/32 glass assembly	£134.00
EVA4132	000432	100mL evaporating flask SJ 29/32	£83.00
EVA1062	000434	500mL evaporating flask SJ 29/32	£83.00
EVA1064	000435	1000mL evaporating flask SJ 29/32	£83.00
EVA1066	20728	1000mL receiving flask SJ 29/32 with P+G safety coating	£137.00
EVA1068	11061167	Foam sensor	£502.00
EVA1070	11059225	AutoDest sensor	£256.00
EVA1072	11060954	Level sensor	£285.00
EVA1074	11061402	Protection shield	£582.00
EVA1076	11059905	Woulff bottle	£196.00
11069167	11069167	Vacuum seal VS26 PTFE base NBR o-ring FDA compliant	£78.00



Buchi V-300 Vacuum Pumps

BUCHI

Chemically resistant PTFE diaphragm vacuum pump with very silent and economical operation. It is optimally suited to be combined with the Rotavapor® R-300.

- High chemical resistance due to the use of high quality materials
- Easy visual detection of condensate and particles due to transparent front panel
- Quiet operation thanks to speed regulation
- Optional post-pump, secondary condenser captures vapour which is otherwise released into the lab atmosphere
- Supplied with vacuum tubing, power cable, silencer, tubing connections, communication cable (Mini-DIN), and membrane tool (to fasten and remove membrane)

Specifications	
Suction Capacity, m ³ /h	1.8
Final Vacuum (absolute), mbar	5 (±2)
Maximum Speed, RPM	1500
Dims, w x d x h, mm	200 x 321 x 291

Code	Alt Ref	Description	Price
EVA1078	11V300100	V-300 with Woulff bottle	£2570.00
EVA1080	11V300010	V-300 with secondary condenser with insulation	£3184.00

A world of products in one place...

- Over 500,000 products
- Easy ordering capability
- Detailed product information
- Promotions and offers

