

## Cling Film

Heavy duty cling film supplied in a dispenser box with cutter.

Code	Dims, mm x m	Price
FIL1003	300 x 300	£31.44
FIL1002	450 x 300	£43.40



## SARAN® Watertight Wrap

Watertight film for the quick closing of containers. No folds or glue. Protects the user from toxic substances or contaminants. 10µm thickness.

Code	Dims, mm x m	Price
FIL1040	300 x 300	£112.00
FIL1042	450 x 300	£92.00



## Parafilm® M Sealing Film

- Protect the integrity of research applications, reduce the risk of tampered data, and keep moisture loss to a minimum
- Odourless, colourless and semi-transparent
- Convenient, flexible and self-sealing, clings around irregular shapes and surfaces
- Stretches more than 200% of original length
- Resistant up to 48 hours against many polar substances such as saline solutions, inorganic acids and alkaline solutions

Code	Dims, mm x m	Price
FIL1020	50 x 75	£55.00
FIL1022	100 x 38	£61.00
FIL1024	100 x 75	£99.00
FIL1026	500 x 15	£119.00



## Heathrow Scientific Parafilm® M Sealing Film Box Top Dispensers



- Accommodates 2 x 50mm rolls or 1 x 100mm rolls
- Easy application, simply place the dispenser on top of the Parafilm® M box and pull the sealing film through, there is no need to remove rolls from the packaging
- Built in safety blade slides and cuts with no danger to fingertips and eliminates the need for scissors
- Supplied with two safety blades, additional replacement blades are available for purchase
- Made from ABS plastic making the unit easy to clean
- Dims, w x d x h: 123 x 172 x 99mm



NEW

Code	Alt Ref	Description	Pack	Price
FIL1270	120667	Blue	Each	£47.00
FIL1272	120758	White	Each	£47.00
FIL1274	120668	Purple	Each	£47.00
FIL1276	120722	Replacement blades	3	£22.00