



SLS Select Folded Filter Discs, Grade Z



Thick creped paper, wet strengthened, used for filtration of fatty products and suitable for hot filtrations. Folded discs.

Supplied in packs of 100

Code	Ø, mm	Price
SLS2133	250	£84.00
SLS2134	330	£100.00



SLS Select Glass Microfibre Filter Discs



Grades GF/1, GF/3 and GF/6 are equivalent to Whatman grades GF/A, GF/C and GF/F respectively.

Supplied in packs of 100

Code	Grade	Pore Size, µm	Ø, mm	Price
SLS4104	GF/1	1.6	25	£17.30
SLS4112	GF/1	1.6	47	£23.80
SLS4120	GF/1	1.6	70	£30.60
SLS4122	GF/1	1.6	90	£45.60
SLS4124	GF/1	1.6	110	£61.00
SLS4126	GF/1	1.6	125	£77.00
SLS4128	GF/1	1.6	150	£95.00
SLS4304	GF/3	1.2	25	£23.40
SLS4310	GF/3	1.2	47	£34.40
SLS4314	GF/3	1.2	70	£39.90
SLS4316	GF/3	1.2	90	£57.00
SLS4318	GF/3	1.2	110	£71.00
SLS4320	GF/3	1.2	125	£104.00
SLS4322	GF/3	1.2	150	£132.00
SLS4514	GF/6	0.7	70	£88.00
SLS4504	GF/6	0.7	25	£38.70
SLS4510	GF/6	0.7	47	£44.70
SLS4516	GF/6	0.7	90	£68.00
SLS4518	GF/6	0.7	110	£97.00
SLS4520	GF/6	0.7	125	£121.00
SLS4522	GF/6	0.7	150	£142.00

Visit our dedicated safety subsite

Bringing you easy to navigate categorisation, extensive choice and live pricing. Visit our safety subsite to view our ever increasing range of safety products.

scientificlabs.co.uk/safety (UK & Northern Ireland)

scientificlabs.ie/safety (Republic of Ireland)

scientificlabs.com/safety (International)



SCIENTIFIC
LABORATORY
SUPPLIES

