

Whatman™ Quantitative Filter Papers Continued

cytiva
Whatman



Ashless, Grade 42

200µm thickness, fast flow, 2.5µm nominal retention.

Supplied in packs of 100

Code	Alt Ref	Ø, mm	Price
FIL2224	1442055	55	£17.00
FIL2225	1442070	70	£17.00
FIL2226	1442090	90	£24.00
FIL2227	1442110	110	£33.00
FIL2228	1442125	125	£39.00
FIL2229	1442150	150	£51.00
FIL2230	1442185	185	£76.00
FIL2231	1442240	240	£119.00



Ashless, Grade 43

220µm thickness, medium to fast flow, 16µm nominal retention.

Supplied in packs of 100

Code	Alt Ref	Ø, mm	Price
FIL2235	1443110	110	£34.00
FIL2236	1443125	125	£40.00
FIL2237	1443150	150	£53.00
FIL2238	1443185	185	£78.00



Ashless, Grade 44

180µm thickness, slow to medium flow, 3µm nominal retention.

Supplied in packs of 100

Code	Alt Ref	Ø, mm	Price
FIL2239	1444070	70	£20.00
FIL2240	1444090	90	£24.00
FIL2241	1444110	110	£33.00
FIL2242	1444125	125	£39.00
FIL2243	1444150	150	£51.00
FIL2244	1444185	185	£76.00



Hardened Low Ash, Grade 589/1

180µm thickness, medium to fast flow, 0.01% typical ash content, 4 - 12µm nominal retention.

Supplied in packs of 100

Code	Alt Ref	Ø, mm	Price
FIL7302	10300009	90	£26.00
FIL7304	10300010	110	£37.00
FIL7306	10300011	125	£43.00
FIL7310	10300014	185	£86.00