

Whatman™ Quantitative Filter Papers Continued

cytiva
Whatman




Hardened Ashless, Grade 540

160µm thickness, medium flow, 8µm nominal retention.

Supplied in packs of 100

Code	Alt Ref	Ø, mm	Price
FIL2301	1540055	55	£18.00
FIL2303	1540090	90	£28.00
FIL2304	1540110	110	£37.00
FIL2305	1540125	125	£44.00
FIL2306	1540150	150	£64.00
FIL2307	1540185	185	£99.00



Hardened Ashless, Grade 541

155µm thickness, fast flow, 22µm nominal retention.

Supplied in packs of 100

Code	Alt Ref	Ø, mm	Price
FIL2311	1541055	55	£20.00
FIL2312	1541070	70	£20.00
FIL2313	1541090	90	£28.00
FIL2314	1541110	110	£36.00
FIL2315	1541125	125	£44.00
FIL2316	1541150	150	£64.00
FIL2317	1541185	185	£97.00
FIL2318	1541240	240	£154.00



Hardened Ashless, Grade 542

150µm thickness, slow flow, 2.7µm nominal retention.

Supplied in packs of 100

Code	Alt Ref	Ø, mm	Price
FIL2322	1542055	55	£20.00
FIL2323	1542070	70	£20.00
FIL2324	1542090	90	£28.00
FIL2325	1542110	110	£38.00
FIL2326	1542125	125	£44.00
FIL2327	1542150	150	£64.00
FIL2328	1542185	185	£91.00