

Whatman™ Qualitative Filter Papers Continued

cytiva
Whatman



Standard, Grade 3

390µm thickness, medium flow thick, 6µm nominal retention.

Supplied in packs of 100

Code	Alt Ref	Ø, mm	Price
FIL2072	1003055	55	£14.00
FIL2073	1003070	70	£14.00
FIL2074	1003090	90	£20.00
FIL2075	1003110	110	£24.00
FIL2076	1003125	125	£28.00
FIL2077	1003150	150	£39.00
FIL2078	1003185	185	£52.00



Standard, Grade 4

205µm thickness, very fast, 25µm nominal retention.

Supplied in packs of 100

Code	Alt Ref	Ø, mm	Price
FIL2086	1004055	55	£10.00
FIL2087	1004070	70	£11.00
FIL2088	1004090	90	£12.00
FIL2089	1004110	110	£15.00
FIL2090	1004125	125	£18.00
FIL2091	1004150	150	£25.00
FIL2092	1004185	185	£35.00
FIL2093	1004240	240	£54.00



Standard, Grade 595

150µm thickness, medium fast thin, 12µm nominal retention.

Supplied in packs of 100

Code	Alt Ref	Ø, mm	Price
FIL7030	10311610	110	£12.00
FIL7031	10311611	125	£15.00
FIL7032	10311612	150	£17.00



Standard, Grade 597

180µm thickness, medium fast, 7µm nominal retention.

Supplied in packs of 100

Code	Alt Ref	Ø, mm	Price
FIL7075	10311808	70	£9.00
FIL7077	10311809	90	£10.00
FIL7079	10311810	110	£12.00
FIL7083	10311812	150	£21.00
FIL7085	10311814	185	£29.00
FIL7086	10311820	240	£47.00