

Sartorius Qualitative Filter Papers

SARTORIUS



Grade 292

The most widely used filter paper for routine applications. Medium retention & flow rate. Qualitative comparable to Whatman Grade 1.

Supplied in packs of 100

Code	Alt Ref	Ø, mm	Price
SAR2026	FT-3-205-042	42	£8.00
SAR2028	FT-3-205-055	55	£10.00
SAR2029	FT-3-205-070	70	£9.00
SAR2031	FT-3-205-090	90	£11.00
SAR2032	FT-3-205-110	110	£11.00
SAR2033	FT-3-205-125	125	£16.00
SAR2034	FT-3-205-150	150	£19.00
SAR2035	FT-3-205-185	185	£29.00
SAR2036	FT-3-205-240	240	£46.00



Grade 6 S/N

Wet strengthened, extra thick with creped surface for removal of coarse particles with a fast flow rate - qualitative 30µm. Alternative to Whatman 113.

Supplied in packs of 50

Code	Alt Ref	Ø, mm	Price
SAR2120	FT-3-314-070	70	£6.00
SAR2121	FT-3-314-090	90	£6.00
SAR2122	FT-3-314-110	110	£8.00
SAR2123	FT-3-314-125	125	£10.00
SAR2124	FT-3-314-150	150	£11.00
SAR2125	FT-3-314-185	185	£17.00
SAR2126	FT-3-314-240	240	£20.00

Sartorius Glass Microfibre Filters

SARTORIUS



Grade MGA

Glass microfibre filter grade for the clarification of buffer and reagent solutions with a particle retention of 1.6µm, weight 54g/m² and thickness 0.23mm. Alternative grade to Whatman GF/A.

Supplied in packs of 100

Code	Alt Ref	Ø, mm	Price
SAR4112	FT-3-1101-047	47	£37.00
SAR4114	FT-3-1101-050	50	£35.00
SAR4116	FT-3-1101-055	55	£35.00
SAR4120	FT-3-1101-070	70	£49.00
SAR4122	FT-3-1101-090	90	£62.00
SAR4124	FT-3-1101-110	110	£82.00
SAR4126	FT-3-1101-125	125	£91.00
SAR4128	FT-3-1101-150	150	£113.00