

OHAUS Pioneer® PX Analytical Balances



Combine essential weighing functionality with excellent performance, at an affordable price and is carefully designed for durability and ease of use.

- High accuracy and repeatability for essential weighing applications
- Cast metal lower housing, sub-pan and stainless steel weighing pan
- Durably constructed for versatile, long term use
- 2nd line display for additional information or guidance, a static removal bar for convenient grounding, and USB connectivity
- Other features include: RS232 interface, integral weigh below hook, removable stainless steel pan, security bracket, software lockout menus, stability indicator, software overload/underload indicators, user selectable environmental settings, user brightness settings, auto dim, auto tare, user selectable operating languages
- Pan dims: Ø 90mm



Code	Alt Ref	Capacity, g	Readability, mg	Calibration	Price
BAL1300	PX124/E	120	0.1	External	£1165.00
BAL1302	PX124	120	0.1	Internal	£1335.00
BAL1304	PX124M	120	0.1	Internal	£1515.00
BAL1306	PX224/E	220	0.1	External	£1280.00
BAL1308	PX224	220	0.1	Internal	£1445.00
BAL1310	PX224M	220	0.1	Internal	£1640.00

Accessories

Code	Alt Ref	Description	Price
BAL1340	80251396	Auxiliary display PAD	£64.00
BAL1342	83034024	Calibrated sinker for liquid density	£305.00
BAL1344	80500525	RS232 cable (9 pin)	£80.70
BAL5146	83021085	USB interface cable	£17.30
BAL5152	80850043	Security device	£157.50
BAL1021	30064202	Printer SF40A	£466.00
BAL5524	80253384	Density kit	£510.00

OHAUS PR Analytical Balances



Designed for routine weighing applications in your workplace offering basic functionality for routine weighing applications. Designed for uncomplicated operation with easy to use display and interface with smart design and durable construction.

- 4 built in applications
- Backlight LCD display
- Constructed with ABS housing, removable stainless steel pan, glass draftshield with sliding top door, integrated weigh-below-hook, security bracket, and calibration lock
- Design features: user selectable environmental filters and brightness settings, auto-tare, autodim, user selectable span calibration points, software lockout and reset menu, user selectable communication settings and data print options, software overload/underload indicator, stability indicator
- Pans dims: Ø 90mm



Code	Alt Ref	Capacity, g	Readability, mg	Calibration	Price
BAL1480	PR124/E	120	0.1	External	£870.00
BAL1482	PR124	120	0.1	Internal	£1015.00
BAL1484	PR224/E	220	0.1	External	£950.00
BAL1486	PR224	220	0.1	Internal	£1095.00